

POWER  
SAFETY  
ACCURACY  
NO COMPROMISE.



## Surgic Pro+

ORAL SURGERY  
MICROMOTOR

**POWERFUL**  
UP TO 80Ncm

**PLUS  
3 FREE  
Accessories**



POWER • SAFETY • ACCURACY

## Surgic Pro+

MODEL	ORDER CODE	RRP PRICE
Surgic Pro Optic	Y1001933	<b>\$8,079</b>

### INCLUDES


- Implant handpiece  
**TiMax X-SG20L**
- Micromotor  
**SGL70M**

**SAVE  
\$2,325**



**FREE  
Travel Case**

ORDER CODE Y1001952  
RRP-\$702



**FREE  
Sterilisation  
Cassette**

MODEL **SG-CASE**  
ORDER CODE Z313101  
RRP-\$584




**FREE  
Irrigation Tube**

BOX OF 50  
ORDER CODE Y900113  
RRP-\$1,039

## Surgical Micromotor & Cable

**FREE  
Sterilisation  
Cassette**



Optic	ORDER CODE	RRP PRICE
SGL70M (with 2m Cord)	E1023001	<b>\$4,133</b>
SGL70M 3.5m (with 3.5m Cord)	E1023002	<b>\$4,498</b>
SG-CASE	Z313101	<b>FREE</b> RRP-\$584

Non-Optic	ORDER CODE	RRP PRICE
SG70M (with 2m Cord)	E1025001	<b>\$3,285</b>
SG70M 3.5m (with 3.5m Cord)	E1025002	<b>\$3,650</b>
SG-CASE	Z313101	<b>FREE</b> RRP-\$584

## iCart Duo

The NSK iCart Duo allows you to install both the NSK SurgicPro and Variosurg Systems as well as additional accessories on to the one easily transportable cart.

MODEL	ORDER CODE	PROMO PRICE
iCart Duo	S9090	<b>\$1,530</b> RRP-\$1,700

### Specifications

- Dimensions : H101.65 cm
- Weight : 16.5 kg
- Adjustable Height
- Lockable Wheels

**10%  
OFF**



"As my business grows NSK has always supported my practice."



My Surgic Pro is consistently reliable, created with precision, and a pleasure to work with.

It is the workhorse of my practice and I couldn't live without it!"

**Dr Cecilia So** BDS Hons (USyd), MFDS (Lond), Grad Dip Clinical Dentistry (Oral Implants) (USyd), FICD, Clinical Associate Lecturer (Oral Implants) (USyd)

## Osseo 100 + Five FREE MultiPegs

### OSSEOINTEGRATION MONITORING DEVICE

MODEL	ORDER CODE	RRP PRICE
Osseo 100	Y1004175	<b>\$3,025</b>




**SAVE  
\$300**

NSK's Osseo 100 measures implant stability and osseointegration to enhance decisions about when to load the implant. This is especially important when using protocols with shorter treatment time and treating higher risk patients.

The uncomplicated operation that measures ISQ, allows the implant loading period to be planned in advance.

## 3-STEP PROCEDURE



The MultiPeg™ is attached to the implant. It screws effortlessly into the implant's internal threads (approximately 6-8 Ncm of torque).



Just aim for the magnet on top of the MultiPeg™. Non-invasive, objective, accurate and repeatable. The peg is excited by magnetic pulses and vibrates due to the stiffness in the contact area between the bone and the implant surface.



An ISQ value is generated and shown on the display. This reflects the level of stability on the universal ISQ scale – from 1 to 99. The higher the ISQ value, the more stable the implant.

## MulTipeg™ REUSABLE

- MulTipeg™
- For all major implant systems \*
- Tissue friendly, durable titanium
- Autoclavable appx. 20 times
- Optimal platform fit
- ISQ Standard Calibrated

\*There are different MulTipegs™ available made to fit different implant system and types. Please refer to the MulTipeg list from the supplier, or scan the QR code to view online.

MODEL	ORDER CODE	PROMO PRICE
MulTipeg™	Various	<b>FREE</b> RRP \$60 EACH



Scan the QR code  
for multipeg list





## Implant Handpieces



MODEL **Ti-Max X-SG20L**  
ORDER CODE C1003

- Optic
- Reducing 20:1
- Titanium
- CA Burs

**PROMO PRICE**  
**\$2,014**

RRP \$2,518



MODEL **Ti-Max X-DSG20L**  
ORDER CODE C1068

- Optic
- Reducing 20:1
- Titanium
- Dismantable
- CA Burs

**PROMO PRICE**  
**\$2,171**

RRP \$2,714



MODEL **S-Max SG20**  
ORDER CODE C1010

- Non Optic
- Reducing 20:1
- Stainless Steel
- CA Burs

**PROMO PRICE**  
**\$1,227**

RRP \$1,534

**20%  
OFF**

## Speed Increasing 1:3



MODEL **Ti-Max X-SG93L**  
ORDER CODE C1004

- Optic
- Increasing 1:3
- Triple spray
- FG Burs

**PROMO PRICE**  
**\$1,935**

RRP \$2,419



MODEL **Ti-Max Z-SG45L**  
ORDER CODE C1107

- Optic
- Increasing 1:3
- Triple Spray
- 25mm FG Burs

**PROMO PRICE**  
**\$1,645**

RRP \$2,056

## Direct Drive 1:1



MODEL **Ti-Max X-SG25L**  
ORDER CODE C1011

- Optic
- Direct Drive 1:1
- CA Burs

**PROMO PRICE**  
**\$1,782**

RRP \$2,228



MODEL **Ti-Max X-SG65L**  
ORDER CODE H1009

- Optic
- Direct Drive 1:1
- CA Burs

**PROMO PRICE**  
**\$1,756**

RRP \$2,195

## Microsurgical Handpieces



MODEL **SGS-ES**  
ORDER CODE H264

- Non Optic
- Direct Drive 1:1
- 65mm HP Burs

**PROMO PRICE**  
**\$1,125**

RRP \$1,407



MODEL **SGA-E2S**  
ORDER CODE H265

- Non Optic
- Increasing 1:2
- 65mm HP Burs

**PROMO PRICE**  
**\$1,830**

RRP \$2,288



MODEL **SGA-ES**  
ORDER CODE H263

- Non Optic
- Direct Drive 1:1
- 20° angulation
- 65mm HP Burs

**PROMO PRICE**  
**\$1,561**

RRP \$1,951



MODEL **SGS-E2S**  
ORDER CODE H266

- Non Optic
- Increasing 1:2
- 20° angulation
- 65mm HP Burs

**PROMO PRICE**  
**\$1,319**

RRP \$1,649

## ISD900

### Faster and Safer Implant Treatments

NSK's iSD900 cordless screwdriver helps to safely place and remove cover screws, healing caps and abutments during dental implant procedures, making treatment up to 50% faster.

MODEL	ORDER CODE	PROMO PRICE
ISD900	Y1001358	<b>\$2,669</b>
		RRP \$3,337

- Universal to major implant systems
- Wide torque range 10-40 Ncm
- Adjustable in 1 or 5 Ncm increments
- Speed range 15, 20, 25 min<sup>-1</sup>



**iSD900**

## DualDefense™ Proprietary Infection Control by NSK

### Clean Head System

Anti-suck back system prevent cross contamination of the air lines

### Anti-Retraction Valve

Anti-Retraction systems prevents cross contamination of the water lines in High-Speed Air Turbines and Electric Micromotors.

## Ti-Max Z **LED** Premium Titanium 2 YEAR WARRANTY

20%  
OFF



MODEL **Ti-Max Z900L**  
ORDER CODE P1111  
• 26W Power  
• Standard Head

PROMO PRICE  
**\$1,467**

RRP \$1,834

MODEL **Ti-Max Z800L**  
ORDER CODE P1110  
• 23W Power  
• Miniature Head

PROMO PRICE  
**\$1,467**

RRP \$1,834

## S-Max M Advanced Stainless Steel 1 YEAR WARRANTY



MODEL **S-Max M900L**  
ORDER CODE P1254  
• 26W Power  
• Standard Head

PROMO PRICE  
**\$1,023**

RRP \$1,279

MODEL **S-Max M800L**  
ORDER CODE P1253  
• 23W Power  
• Miniature Head

PROMO PRICE  
**\$1,023**

RRP \$1,279



## S-Max Pico Paediatric Highspeed

ULTRA MINI HEAD  
1 YEAR WARRANTY

MODEL **S-Max Pico**  
ORDER CODE P1140

- 8W Power
- Ultra Short Burs available

PROMO PRICE  
**\$1,023**

RRP \$1,279

Non-Optic options available.  
Other connections offered for KaVo, Sirona, W&H, Bien-Air.



## Ti-Max X450L Surgical Highspeed

REAR EXHAUSTED AIR  
2 YEAR WARRANTY

MODEL **Ti-Max X450L**  
ORDER CODE P1078

- 21W Power
- 45° head
- Triple Water Jets

PROMO PRICE  
**\$1,242**

RRP \$1,562

## Buy ANY 6 Ti Max Handpieces Get **FREE** iCare\*

\*Buy any 6 Ti Max handpieces at RRP and receive a FREE iCare and adaptors of your choice!

- Pre-set solution volume based on handpiece type
- Thorough cleaning and lubrication with simple operation
- Chuck cleaning and lubricating function
- Purging mode removes excess oil
- Integrated 1.2 litre oil tank reduces burden of refill

### iCare C2 Type

ORDER CODE Y1002787

Dual high speed + Dual low speed  
Connectors requires 2 high  
speed Adaptors\*



### iCare C3 Type

ORDER CODE Y1002788

Single high speed + Triple low  
speed Connectors requires  
1 high speed Adaptor\*



\*Adaptors available to fit NSK, KaVo, Sirona, Bien-Air, Morita high speed handpieces.  
Limit x1 free adaptor free per each iCare C3 and x2 free adaptors per each iCare C2.



**iCare**



## iClave plus



**FREE**  
BAR CODE PRINTER  
& WI-FI LOG  
SAVE \$1,569



Bar Code Printer



Wi-Fi Log

NEW

\* Available via authorised NSK Autoclave dealers

## iClave Plus + **FREE** accessories

B-type vacuum autoclave\*

MODEL	ORDER CODE	RRP PRICE
iClave plus	Y1003077	<b>\$9,990</b>

- 20% more sterilisation space than a conventional stainless steel chamber of the same size.
- Copper chamber minimises internal temperature fluctuations.
- Super fast cycle as quick as 18 minutes, including drying.
- One year warranty.
- Unit size: H428 x W445 x D532mm, chamber 18L.

## Optional PURITY Water Purifier

MODEL	ORDER CODE	PROMO PRICE
PURITY OR2	Z1284001	<b>\$1,301</b>

Plus free shipping if buying both offers!

RRP \$2,602



**50%  
OFF**

**NSK Oceania Pty Ltd** [www.nskoceania.com.au](http://www.nskoceania.com.au)

Product Information | Technical Support

Freephone: 1300 44 33 21 (Australia) Email: [info@nskoceania.com.au](mailto:info@nskoceania.com.au)