

# NSK is going **green**

Renew. Recycle. Reuse.



# HANDPIECES AND MICROMOTORS

## High Speed Handpieces

### S-Max and Ti-Max Series



MODEL **Ti-Max Z800L**  
ORDER CODE **P1110**

- Cellular Glass Optics
- Quattro Spray
- Miniature Head

**TRADE-IN PRICE**  
**\$1,376**  
**SALE PRICE**  
**\$1,559**  
RRP \$1,834



MODEL **Ti-Max Z900L**  
ORDER CODE **P1111**

- Cellular Glass Optics
- Quattro Spray
- Miniature Head

**TRADE-IN PRICE**  
**\$1,376**  
**SALE PRICE**  
**\$1,559**  
RRP \$1,834



MODEL **S-Max M800L**  
ORDER CODE **P1253**

- 23W Power
- Quattro Spray
- Miniature Head

**TRADE-IN PRICE**  
**\$959**  
**SALE PRICE**  
**\$1,087**  
RRP \$1,279



MODEL **S-Max M900L**  
ORDER CODE **P1254**

- 26W Power
- Quattro Spray
- Standard Head

**TRADE-IN PRICE**  
**\$959**  
**SALE PRICE**  
**\$1,087**  
RRP \$1,279



MODEL **Ti-Max Z95LSW**  
ORDER CODE **C1034SW**

- Optic
- 1:5 Increasing
- Quattro Spray
- Titanium Body Coating DURAGRIP

**TRADE-IN PRICE**  
**\$1,473**  
**SALE PRICE**  
**\$1,669**  
RRP \$1,964



MODEL **Ti-Max Z25L**  
ORDER CODE **C1038**

- Optic
- 1:1 Direct Drive
- Single Spray
- For CA burs (ø2.35)

**TRADE-IN PRICE**  
**\$1,035**  
**SALE PRICE**  
**\$1,173**  
RRP \$1,380



MODEL **S-Max M95L**  
ORDER CODE **C1023**

- Optic
- 1:5 Increasing
- Quattro Spray
- For FG burs (ø1.6)

**TRADE-IN PRICE**  
**\$1,224**  
**SALE PRICE**  
**\$1,387**  
RRP \$1,632



MODEL **S-Max M25L**  
ORDER CODE **C1024**

- Optic
- 1:1 Direct Drive
- Single Spray
- For CA burs (ø2.35)

**TRADE-IN PRICE**  
**\$939**  
**SALE PRICE**  
**\$1,064**  
RRP \$1,252



MODEL **S-Max Pico**  
ORDER CODE **P1140**

- Optic
- Ultra Mini Head
- 8W Power
- Ultra Short Burs available

**TRADE-IN PRICE**  
**\$959**  
**SALE PRICE**  
**\$1,087**  
RRP \$1,279



MODEL **Ti-Max X450L**  
ORDER CODE **P1078**

- Optic
- Rear Exhausted Air
- 21W Power
- 45° Head
- Triple Spray

**TRADE-IN PRICE**  
**\$1,164**  
**SALE PRICE**  
**\$1,319**  
RRP \$1,552



MODEL **Ti-Max X65L**  
ORDER CODE **H273**

- 1:1 Direct drive
- Clean Head System
- Single spray
- For HP burs (ø2.35)
- For CA burs (ø2.35)

**TRADE-IN PRICE**  
**\$1,357**  
**SALE PRICE**  
**\$1,538**  
RRP \$1,810



MODEL **S-Max M65**  
ORDER CODE **H1008**

- Non-Optic
- 1:1 Direct Drive
- Single Spray
- HP/CA burs

**TRADE-IN PRICE**  
**\$567**  
**SALE PRICE**  
**\$643**  
RRP \$757



MODEL **S-Max M205**  
ORDER CODE **M1007001**

- Non-Optic
- 1:1 Direct Drive
- Internal Spray
- Non retraction valve

**TRADE-IN PRICE**  
**\$843**  
**SALE PRICE**  
**\$955**  
RRP \$1,124



MODEL **Ti-Max X205L**  
ORDER CODE **M1013**

- Optic LED
- 1:1 Direct Drive
- Internal Spray
- 360° Swivel
- 32,000 LUX

**TRADE-IN PRICE**  
**\$1,365**  
**SALE PRICE**  
**\$1,547**  
RRP \$1,820

Non-optic options available.  
Other connections offered for KaVo, Sirona, W&H, Bien Air.

## Going green changes lives

NSK's Going Green campaign collects, refurbishes and repairs used and non-performing handpieces and equipment which are destined for landfill. The handpieces have so far been distributed to areas including charities like "Smiles for the Pacific" which support various Polynesian Island dental hospitals.

Photo: NSK Oceania Managing Director Phil Fitzgibbon (third right) unites with members of the Solomon Island Oral Health Clinic as part of the Smiles for the Pacific program – one of NSK's initiatives to promote access to dentistry.



## FX Series (NON-OPTIC)



MODEL **FX15M**  
ORDER CODE **C1053001**

- 4:1 Reduction
- Latch type prophyl cups
- For CA burs (ø2.35)
- Max Speed: 10,000 min<sup>-1</sup>

**TRADE-IN PRICE**  
**\$565**  
**SALE PRICE**  
**\$597**  
RRP \$628



MODEL **FX25M**  
ORDER CODE **C1052001**

- 1:1 Direct drive
- For CA burs (ø2.35)
- Push button chuck
- Ball bearing
- Max Speed: 40,000 min<sup>-1</sup>

**TRADE-IN PRICE**  
**\$385**  
**SALE PRICE**  
**\$406**  
RRP \$427



MODEL **FX57M**  
ORDER CODE **C1054001**

- 4:1 Reduction
- For screw-in cups and brushes
- 360° Rotation
- Paste Defense System
- Max Speed: 5,000 min<sup>-1</sup>

**TRADE-IN PRICE**  
**\$423**  
**SALE PRICE**  
**\$447**  
RRP \$470



MODEL **FX205M M4**  
ORDER CODE **M1005001**

- 1:1 Direct Drive
- External Spray
- Enhanced Balance
- Powerful Torque

**TRADE-IN PRICE**  
**\$884**  
**SALE PRICE**  
**\$933**  
RRP \$982



This offer is for one used item of any brand, for one equivalent new NSK item. Not valid for parts. Items must be whole but not necessarily in working condition. Valid until 31st December 2021. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.

## Surgic Pro

### Implant and Oral Surgery

Multi-Functional Brushless Optic Surgical Micromotor for Oral Surgery and Implantology.



- LED Illumination
- Powerful torque (up to 80Ncm)
- Wide speed range (200 to 40,000 min<sup>-1</sup>)
- High torque accuracy
- Lightest and smallest micromotor
- LED illumination (over 32,000 LUX)
- Ti Max X-SG20L implant handpiece

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
Surgic Pro	Y1001933	\$6,059	\$6,867
			RRP \$8,079.23

## Osseo 100

Osseo 100 measures implant stability and osseointegration using RFA, to enhance decisions about when to load the implant.

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
Osseo 100	Y1004175	\$2,269	\$2,571
			RRP \$3,025



**25% OFF**  
WITH TRADE-IN PROMO

**15% OFF**  
WITH OCT-DEC PROMO

## ISD900

### Faster and Safer Implant Treatments

NSK's iSD900 cordless screwdriver helps to safely place and remove cover screws, healing caps and abutments during dental implant procedures, making treatment up to 50% faster.

- Universal to major implant systems
- Wide torque range 10-40 Ncm
- Adjustable in 1 or 5 Ncm increments
- Speed range 15, 20, 25 min<sup>-1</sup>

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
ISD90	Y1001358	\$2,275	\$2,578
			RRP \$3,033



## Surgical Handpieces

### Implant Handpieces



MODEL **Ti-Max X-SG20L**  
ORDER CODE C1003

- Optic
- Reducing 20:1
- Titanium
- CA Burs

TRADE-IN PRICE  
**\$1,888**

SALE PRICE  
**\$2,140**

RRP \$2,518



MODEL **S-Max SG20**  
ORDER CODE C1010

- Non Optic
- Reducing 20:1
- Stainless Steel
- CA Burs

TRADE-IN PRICE  
**\$1,150**

SALE PRICE  
**\$1,303**

RRP \$1,534

### Speed Increasing 1:3



MODEL **Ti-Max X-SG93L**  
ORDER CODE C1004

- Optic
- Increasing 1:3
- Triple spray
- FG Burs

TRADE-IN PRICE  
**\$1,814**

SALE PRICE  
**\$2,056**

RRP \$2,419



MODEL **Ti-Max Z-SG45L**  
ORDER CODE C1107

- Optic
- Increasing 1:3
- Triple Spray
- 25mm FG Burs

TRADE-IN PRICE  
**\$1,542**

SALE PRICE  
**\$1,748**

RRP \$2,056

### Direct Drive 1:1



MODEL **Ti-Max X-SG25L**  
ORDER CODE C1011

- Optic
- Direct Drive 1:1
- CA Burs

TRADE-IN PRICE  
**\$1,671**

SALE PRICE  
**\$1,894**

RRP \$2,228



MODEL **Ti-Max X-SG65L**  
ORDER CODE H1009

- Optic
- Direct Drive 1:1
- CA Burs

TRADE-IN PRICE  
**\$1,646**

SALE PRICE  
**\$1,865**

RRP \$2,195

## Microsurgical Handpieces



MODEL **SGS-ES**  
ORDER CODE H264

- Non Optic
- Direct Drive 1:1
- 65mm HP Burs

TRADE-IN PRICE  
**\$1,055**

SALE PRICE  
**\$1,196**

RRP \$1,407



MODEL **SGA-E2S**  
ORDER CODE H265

- Non Optic
- Increasing 1:2
- 65mm HP Burs

TRADE-IN PRICE  
**\$1,716**

SALE PRICE  
**\$1,945**

RRP \$2,288



MODEL **SGA-ES**  
ORDER CODE H263

- Non Optic
- Direct Drive 1:1
- 20° angulation
- 65mm HP Burs

TRADE-IN PRICE  
**\$1,463**

SALE PRICE  
**\$1,658**

RRP \$1,951



MODEL **SGS-E2S**  
ORDER CODE H266

- Non Optic
- Increasing 1:2
- 20° angulation
- 65mm HP Burs

TRADE-IN PRICE  
**\$1,237**

SALE PRICE  
**\$1,402**

RRP \$1,649

This offer is for one used item of any brand, for one equivalent new NSK item. Not valid for parts. Items must be whole but not necessarily in working condition. Valid until 31st December 2021. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.



## Varios Combi Pro Basic Kit

### 2 Maintenance Systems in 1

Varios Combi Pro combines ultrasonic scaling and powder therapy for a more effective and comfortable oral hygiene and biofilm control.

#### Basic Kit

MODEL **Varios Combi Pro Basic Kit**  
ORDER CODE **Y1004194**

- Advanced optic iPiezo scaling for General, Perio and Endo treatment.
- Powder Therapy to clean Supragingival areas efficiently and quickly.
- User friendly handpiece and nozzle for easy access.

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
Varios Combi Pro Basic Kit	Y1004194	<b>\$5,455</b>	<b>\$6,183</b> RRP \$7,274



Scan the QR code  
for more information



## Varios 970

### Portable/Bench Top

Optimum clinically versatile self-contained portable bench top unit, able to operate using two 400ml solution bottles with independent pumps and an optional water connection.

- iPiezo engine ensures low vibration and stabilised power
- Comfortable, lightweight handpiece
- Easy switch from Endo to Perio to Scaling and General
- LED illumination
- Memory function and LCD display

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
Optic Complete Set	Y1001168	<b>\$3,143</b>	<b>\$3,562</b> RRP \$4,190



## Prophy-Mate neo

### Tooth Polishing System

Prophy-Mate neo, the best polishing tool for beautiful white teeth.

- Removal of plaque and discolouration
- Etching of enamel, dentine and porcelain prior to bonding
- Quick connect to major brands of turbine couplings
- Preparation for pit and fissure sealant
- Sturdy and reliable
- Simple to operate
- Easy maintenance
- Available with 60° or 80° nozzles




MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
PMNG-PTL-P for NSK Coupling	Y135029	<b>\$1,177</b>	<b>\$1,334</b> RRP \$1,570

## Perio-Mate

### Biofilm Eraser

Perio-Mate, the best biofilm eraser for gentle periodontal treatment.

- Treatment of periodontal pockets and peri-implantitis
- Eradication of biofilm in periodontal pockets
- Removal of plaque and discolouration
- NSK's innovative nozzle tip
- Precise, flawless powder flow
- Slim nozzle design
- Smooth 360° rotation
- Powder Control Ring



MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
PER-PTL-P for NSK Coupling	Y1002654	<b>\$1,177</b>	<b>\$1,334</b> RRP \$1,570

This offer is for one used item of any brand, for one equivalent new NSK item. Not valid for parts. Items must be whole but not necessarily in working condition.  
Valid until 31st December 2021. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.

## LABORATORY

### Volvere i7

#### Coreless Micromotor

- Light and compact to fit into the smallest space
- Speed 1,000-35,000 min<sup>-1</sup>
- Microcomputer Control
- Auto Cruise Function
- Contents: Control unit, Micromotor Handpiece, Foot control, Handpiece Stand

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
Volvere i7 RM	Y1003493	<b>\$1,468</b>	<b>\$1,664</b>
			RRP \$1,958



Standard Micromotor



E-Type Micromotor

### Ultimate XL

#### Bench Top System

- Speed 1,000-50,000 min<sup>-1</sup>
- Max Torque 8.7Ncm
- Precise control by microcomputer
- User-friendly design
- Auto Cruise Function
- Easy-to-operate Display Panel

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
Ultimate XL Torque micromotor: UMXL-GT (230V)	Y141469	<b>\$2,036</b>	<b>\$2,307</b>
			RRP \$2,714



## ENDODONTICS

### Simple and precise root canal treatment

#### Endodontic Motors

##### Endo-Mate AT

MODEL Endo-Mate AT (230V)  
ORDER CODE Y1002190

- Fits all major brands of files
- Torque control, with integrated auto-reverse function
- Compact stylish control/monitor pad for simple and user friendly operation

TRADE-IN PRICE  
**\$1,256**  
SALE PRICE  
**\$1,424**

RRP \$1,675



##### Endo-Mate TC2

MODEL Endo-Mate TC2 (230V)  
MP-F16R Head

ORDER CODE Y1001027

- Large LCD display with 5 programs for different file systems
- Supports most major brands of Ni-Ti files
- Torque control facility ranging from 0.1-4.0 NCM depending on head:gear ratio

TRADE-IN PRICE  
**\$1,571**  
SALE PRICE  
**\$1,781**

RRP \$2,095



#### Apex Locator

##### iPex II

MODEL iPEX II Complete Set  
ORDER CODE Y1002208

- Accurate measurement of any root canal
- SmartLogic technology
- Operates in all conditions, including dry, wet and bloody canals

TRADE-IN PRICE  
**\$1,177**  
SALE PRICE  
**\$1,334**

RRP \$1,570



This offer is for one used item of any brand, for one equivalent new NSK item. Not valid for parts. Items must be whole but not necessarily in working condition. Valid until 31st December 2021. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.

## iCare C2/C3

### Handpiece oiling and maintenance system

- Pre-set solution volume based on handpiece type
- Thorough cleaning and lubrication with simple operation
- Chuck cleaning and lubricating function
- Purging mode removes excess oil
- Integrated 1.2 litre oil tank reduces burden of refill



MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
<b>iCare C2</b> Standard System Complete • 2 rotation ports • 2 non-rotation ports • 230V	Y1002787	<b>\$3,812</b>	<b>\$4,320</b> <small>RRP \$5,082</small>
<b>iCare C3</b> Standard System Complete • 3 rotation ports • 1 non-rotation port • 230V	Y1002788	<b>\$4,071</b>	<b>\$4,614</b> <small>RRP \$5,429</small>



## Adaptors

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
Adaptor for NSK High Speed Handpieces that fit NSK FlexiQuik type coupling	Z257010	<b>\$165</b>	<b>\$188</b> <small>RRP \$221</small>
Adaptor for Kavo High and Low Speed that fit Kavo Multiflex Coupling	Z257020	<b>\$165</b>	<b>\$188</b> <small>RRP \$221</small>
Adaptor for W&H High Speed to fit W&H Roto Quick Coupling	Z257021	<b>\$165</b>	<b>\$188</b> <small>RRP \$221</small>
Adaptor for Bien Air High and Low Speed that fit Bien Air Unifix Coupling	Z257022	<b>\$165</b>	<b>\$188</b> <small>RRP \$221</small>
Adaptor for Sirona High and Low Speed that fit Sirona Quick Coupling	Z257023	<b>\$165</b>	<b>\$188</b> <small>RRP \$221</small>
Adaptor For Endo-Mate DT/TC, Taskal Wizards Heads	Z257042	<b>\$165</b>	<b>\$188</b> <small>RRP \$221</small>
Adaptor to fit Morita High Speed Handpiece	Z257048	<b>\$165</b>	<b>\$188</b> <small>RRP \$221</small>



## iClave plus

### B-type vacuum autoclave

- 18 litre chamber capacity.
- 20% more sterilisation space than a conventional stainless steel chamber of the same size.
- Copper chamber minimises internal temperature fluctuations.
- Super fast cycle as quick as 18 minutes, including drying.
- One year warranty.

Unit size: H428 x W445 x D532mm, chamber 18L.

MODEL	ORDER CODE	TRADE-IN PRICE	SALE PRICE
<b>iClave plus</b>	Y1003077	<b>\$8,465</b>	<b>\$9,405</b> <small>RRP \$9,900</small>



**NSK Oceania Pty Ltd** [www.nskoceania.com.au](http://www.nskoceania.com.au)

Product Information | Technical Support

Freephone: 1300 44 33 21 (Australia) Email: [info@nskoceania.com.au](mailto:info@nskoceania.com.au)

This offer is for one used item of any brand, for one equivalent new NSK item. Not valid for parts. Items must be whole but not necessarily in working condition. Valid until 31st December 2021. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.