

Meet nano.

NEW

Ti-Max nano

Contra-Angle

Small, but mighty.


Delicate and precise, a natural extension of your hand.

10% smaller and lighter*



Manufacturer Warranty

Means peace of mind.



NSK Products are precision manufactured in Japan with 100% of NSK parts manufactured in-house. It is with this quality in mind that NSK products come with guarantees of the highest standard. Registering your NSK products today for Warranty will ensure total peace of mind. With NSK, you are supported for the life of your NSK products.



To register your product for warranty visit **NSKtech.com.au** or simply scan the QR Code

Your local NSK Business Development and Education Specialists

WA
Stuart Blair
0437 407 668
stuart@nskoceania.com.au

VIC, SA & NT
Mariam Dalzotto
0427 833 684
mariam@nskoceania.com.au

Queensland
Sam Vaizie
0439 260 006
sam@nskoceania.com.au

NSW North
Theresa Sparkes
0419 251 216
theresa@nskoceania.com.au

NSW South and ACT
Janette Carroll
0439 199 979
janette@nskoceania.com.au

Tasmania
Stuart Blair
0437 407 668
stuart@nskoceania.com.au

CONTACT US

NSK Oceania Pty Ltd

www.nskoceania.com.au

Unit 12, 809-821 Botany Road, Rosebery NSW 2018, Australia

Product Information | Technical Support
Freephone: 1300 44 33 21



Ti-Max nano


Shorter, Perfectly Balanced

Lighter and perfectly balanced, "Ti-Max nano" handpieces are designed to be a natural extension of the clinician's hand, ensuring smoother procedures. Even during procedures that require delicate and precise movements, a "nano" handpiece faithfully mirrors the user's intent. NSK continues to challenge the limits of manufacturing to pioneer next-generation contra-angles.



Smart Miniaturization Enables Perfect Balance

NSK's Ti-Max nano preserves the high performance that NSK contra-angles are known for. These new handpieces are approximately 10 % smaller and lighter than conventional models and in combination with NLZ micromotor series, together they give you the ease of manoeuvrability and accuracy you are used to from your air-turbine.



Perfect Balance

Reducing size and weight shifts the centre of gravity to an optimal position, dramatically reducing fatigue during procedures and improving your fingertip dexterity.



Clean Head System

The Clean Head System has a labyrinth seal mechanism that reduces blood and foreign matter from entering the handpiece head from the oral cavity. As a result, the ball bearings and mechanical parts last longer.



Smaller and Slimmer

nano95LS has a smaller head and slimmer neck dimensions, which increases accessibility and visibility to posterior and difficult access areas.

Ti-Max nano line up

Exceptionally small head and slim body / Clear visibility and easy access / Amazingly quiet operation / Cellular Glass Optics / Titanium body with scratch resistant DURAGRIP / Clean Head System / Push Button Chuck (exclude nano65LS) / Outstanding Durability.



Optic		
MODEL	ORDER CODE	RRP
nano95LS	C1099	\$1,909
		Sale Price
		\$1,623

Specifications

- Gear Ratio: 1:5 Increasing
- For: FG Burs (Ø1.6)
- Water Spray: Quattro Spray
- Body Material: Titanium
- Body Coating: DURAGRIP
- Max Speed: 200,000 min⁻¹

Features:

- Anti Heat System
- Ceramic Bearings
- Microfilter

Optic		
MODEL	ORDER CODE	RRP
nano15LS	C1101	\$1,907
		Sale Price
		\$1,621

Specifications

- Gear Ratio: 4:1 Reduction
- For: CA Burs (Ø2.35)
- Water Spray: Single Spray
- Body Material: Titanium
- Body Coating: DURAGRIP
- Max Speed: 10,000 min⁻¹

Features:

- Microfilter

Optic		
MODEL	ORDER CODE	RRP
nano25LS	C1100	\$1,342
		Sale Price
		\$1,141

Specifications

- Gear Ratio: 1:1 Direct Drive
- For: CA Burs (Ø2.35)
- Water Spray: Single Spray
- Body Material: Titanium
- Body Coating: DURAGRIP
- Max Speed: 40,000 min⁻¹

Features:


- Microfilter

Optic		
MODEL	ORDER CODE	RRP
nano65LS	H1048	\$1,760
		Sale Price
		\$1,496

Specifications

- Gear Ratio: 1:1 Direct Drive
- For: HP / CA (Ø2.35)*
- Water Spray: Single Spray
- Body Material: Titanium
- Body Coating: DURAGRIP
- Max Speed: 40,000 min⁻¹

*Need to insert the supplied bur stopper



Buy iClave and get 50% off iCare



iClave plus
B-type vacuum autoclave

- 18 litre chamber capacity.
 - 20% more sterilisation space than a conventional stainless steel chamber of the same size.
 - Copper chamber minimises internal temperature fluctuations.
 - Super fast cycle as quick as 18 minutes, including drying.
 - One year warranty.
- Unit size: H428 x W445 x D532mm, chamber 18L.

MODEL	ORDER CODE	RRP
iClave plus	Y1003077	\$10,692

iCare C2/C3
Handpiece oiling and maintenance system

- Pre-set solution volume based on handpiece type.
 - Thorough cleaning and lubrication with simple operation.
 - Chuck cleaning and lubricating function.
 - Purging mode removes excess oil.
 - Integrated 1.2 litre oil tank reduces burden of refill.
- Adaptors sold separately. Available to suit NSK, KaVo, W&H, Sirona, Bien Air & Morita connections.

Model	ORDER CODE	RRP	Buy iClave and GET 50% OFF
iCare C2 Standard System Complete <ul style="list-style-type: none">• 2 rotation ports• 2 non-rotation ports• 230V	Y1002787	\$5,489	\$2,745
iCare C3 Standard System Complete <ul style="list-style-type: none">• 3 rotation ports• 1 non-rotation port• 230V	Y1002788	\$5,863	\$2,932

Set up your surgery with NSK and save

Huge offers on products and packages

T&Cs apply, to find out more contact your local dealer today!

SAVE \$1,000s
On surgery setup costs

REDUCE AEROSOLS GO ELECTRIC

Anti-Retraction
Valve prevents
suckback/
contamination of
the water line



Accurate and clean margins are the most important factor when achieving the best digital scan. Electric driven speed increasing handpieces are the best options when looking for the control you need to produce accurate margins.

What is a speed increasing handpiece?

A contra-angle that does everything a turbine (high speed) does and uses the same FG burs. But it has incredible 'torque' originating from the micromotor which puts it in a league of its own and puts you in total control – not more loss of power as you come into contact with the tooth. Choose your speed. Save time cutting preps.

What motor do I need to use speed increasing handpieces?

It must be an electric micromotor (not an air motor) as the gears inside (usually 1:5) need the electrical power to run at a high speed. These ultimate brushless electric micromotors can be retro fitted by an engineer or come as an option when buying new chairs. NSK's NLX micromotors offer a powerful 3.4Ncm of torque and the new NLZ micromotor offers increased torque up to 4.2 Ncm.



Ti-Max Z

BUY 3, GET 1 FREE*

S-Max M

NSK Ti-Max Z range

	MODEL		ORDER CODE	RRP
Optic	Ti-Max Z95L	1:5 Increasing	C1034SW	\$2,121
Optic	Ti-Max Z85L	1:5 Increasing	C1062	\$2,121
Optic	Ti-Max Z45L	1:4.2 Increasing	C1064	\$2,263
Optic	Ti-Max Z25L	1:1 Direct Drive	C1038	\$1,491
Optic	Ti-Max Z15L	4:1 Reducing	C1039	\$2,119

NSK S-Max M range

	MODEL		ORDER CODE	RRP
Optic	S-Max M95L	1:5 Increasing	C1023	\$1,763
Non-Optic	S-Max M95		C1026	\$1,405
Optic	S-Max M25L	1:1 Direct Drive	C1024	\$1,353
Non-Optic	S-Max M25		C1027	\$949
Optic	S-Max M15L	4:1 Reducing	C1025	\$1,647
Non-Optic	S-Max M15		C1028	\$1,334
Non-Optic	S-Max M65	1:1 Direct Drive	H1008	\$818

Buy Surgic Pro complete and get 50% off VarioSurg3

Surgic Pro

Implant and Oral Surgery

Multi-Functional Brushless Optic Surgical Micromotor for Oral Surgery and Implantology.



- LED Illumination
- Powerful torque (up to 80Ncm)
- Wide speed range (200 to 40,000 min⁻¹)
- High torque accuracy
- Lightest and smallest micromotor
- LED illumination (over 32,000 LUX)
- Ti Max X-SG20L implant handpiece
- Excellent durability
- Autoclavable, thermo-disinfectable and maintenance free brushless design

**BUNDLE
AND SAVE
50%**

MODEL	ORDER CODE	RRP
Surgic Pro	Y1001933	\$8,726

VarioSurg3

Ultrasonic Bone Surgery System

Safer Bone Surgery

VarioSurg3 Ultrasonic bone surgery system provides safe and precise cutting of bone, without damaging soft tissue. A variety of tips makes procedures easier and stress-free to help provide more precise and efficient treatment.



- Major Power Boost
- 3 Modes (SURG, ENDO, PERIO)
- Selectable Burst Mode
- Adjustable LED Light intensity
- Automatic Cleaning Mode
- 9 Programs (SURGx5, ENDOx2, PERIOx2)
- Feedback Function
- Auto Tuning Function
- Dynamic Link Feature

MODEL	ORDER CODE	RRP	Buy Surgic Pro and GET 50% OFF
VarioSurg3	Y1002248	\$10,060	\$5,030

Create Synergy in Implantology by linking the Surgic Pro & Variosurg together.

Linking options

SG Link Set

MODEL	ORDER CODE	RRP
VarioSurg3 SG Link Set (Stand & Cable)	Y1002729	\$1,116



iCart Duo & Link Cable

MODEL	ORDER CODE	RRP
iCart Duo	S9090	\$2,355
Link Cable	Z1189	\$163



Buy Surgic Pro complete and get FREE extra surgical handpiece

**THE
SURGICAL
HANDPIECE
OPTIONS
AVAILABLE**



MODEL **Ti-Max X-SG20L**
ORDER CODE C1003

- Optic
- Reducing 20:1
- Titanium
- CA Burs

RRP
\$2,720

SALE PRICE

FREE



MODEL **Ti-Max Z-SG45L**
ORDER CODE C1107

- Optic
- Increasing 1:3
- Triple Spray
- 25mm FG Burs

RRP
\$2,221

SALE PRICE

FREE



MODEL **Ti-Max X-SG65L**
ORDER CODE H1009

- Optic
- Direct Drive 1:1
- HP Burs

RRP
\$2,371

SALE PRICE

FREE

NSK Oceania Pty Ltd www.nskoceania.com.au

Product Information | Technical Support

Freephone: 1300 44 33 21 (Australia) Email: info@nskoceania.com.au

This offer is for one used item of any brand, for one equivalent new NSK item. Not valid for parts. Items must be whole but not necessarily in working condition. Valid until 31st March 2022. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.