

Full Powered Flexibility

Powerful 26W
The smooth and comfortable usability



Ti-Max Z
Air Turbine

High Speed Handpieces

S-Max and Ti-Max Series



MODEL **Ti-Max Z800L**
ORDER CODE P1110

- Optic
- 23W Power
- Miniature head

RRP	BUY 1
\$1,981	\$1,684
BUY 2	BUY 3
\$1,585	\$1,486



MODEL **Ti-Max Z900L**
ORDER CODE P1111

- Optic
- 26W Power
- Standard head

RRP	BUY 1
\$1,981	\$1,684
BUY 2	BUY 3
\$1,585	\$1,486



MODEL **S-Max M800L**
ORDER CODE P1253

- Optic
- 23W Power
- Miniature head

RRP	BUY 1
\$1,381	\$1,174
BUY 2	BUY 3
\$1,105	\$1,036



MODEL **S-Max M900L**
ORDER CODE P1254

- Optic
- 26W Power
- Standard head

RRP	BUY 1
\$1,381	\$1,174
BUY 2	BUY 3
\$1,105	\$1,036



MODEL **Ti-Max X450L**
ORDER CODE P1078

- Optic
- Re exhausted 45° head
- 21W Power

RRP	BUY 1
\$1,677	\$1,425
BUY 2	BUY 3
\$1,342	\$1,258



MODEL **S-Max Pico**
ORDER CODE P1140

- Optic
- 8W Power
- Ultra Mini Head

RRP	BUY 1
\$1,381	\$1,174
BUY 2	BUY 3
\$1,105	\$1,036



Contra Angle Handpieces

S-Max and Ti-Max Series



MODEL **Ti-Max Z95L**
ORDER CODE C1034SW

- Optic
- 1:5 Speed Increasing
- Switch Spray function
- For FG Burs

RRP	BUY 1
\$2,121	\$1,803
BUY 2	BUY 3
\$1,697	\$1,591



MODEL **Ti-Max Z45L**
ORDER CODE C1064

- Optic
- 1:4.2 Speed Increasing
- Switch Spray function
- 45° Head
- For FG burs

RRP	BUY 1
\$2,263	\$1,924
BUY 2	BUY 3
\$1,810	\$1,697



MODEL **Ti-Max Z85L**
ORDER CODE C1062

- Optic
- 1:5 Speed Increasing
- Single Spray
- For FG Burs

RRP	BUY 1
\$2,221	\$1,803
BUY 2	BUY 3
\$1,697	\$1,591



MODEL **Ti-Max Z25L**
ORDER CODE C1038

- Optic
- 1:1 Direct Drive
- Single spray
- For CA burs (ø2.35)

RRP	BUY 1
\$1,491	\$1,267
BUY 2	BUY 3
\$1,193	\$1,118



MODEL **Ti-Max Z24L**
ORDER CODE C1126001

- Optic
- 1:1 Direct Drive
- Single spray
- For CA burs (ø2.35)

RRP	BUY 1
\$1,491	\$1,267
BUY 2	BUY 3
\$1,193	\$1,118



MODEL **S-Max M95L**
ORDER CODE C1023

- Optic
- 1:5 Increasing
- Quattro Spray
- For FG burs (ø1.6)

RRP	BUY 1
\$1,763	\$1,499
BUY 2	BUY 3
\$1,410	\$1,322



MODEL **S-Max M25L**
ORDER CODE C1024

- Optic
- 1:1 Direct Drive
- Single Spray
- For CA burs

RRP	BUY 1
\$1,353	\$1,150
BUY 2	BUY 3
\$1,082	\$1,015



MODEL **S-Max M15L**
ORDER CODE C1025

- Optic
- 4:1 Reducing
- Single Spray
- For CA burs

RRP	BUY 1
\$1,647	\$1,400
BUY 2	BUY 3
\$1,318	\$1,235

Non-optic options available. Other connections offered for KaVo, Sirona, W&H, Bien Air.

Guaranteed protection. Certified excellence.

The industry leader in infection control for 30 years.

Our proprietary, best-in-class Clean Head System prevents contaminants from entering the handpiece, while our advanced Non-Retraction Valve protects against cross-contamination. NSK sets the standard for safety and protection with the world's most technologically advanced, award-winning handpiece design.



MODEL **Ti-Max X65L**
ORDER CODE H273

- 1:1 Direct drive
- Clean Head System
- Single spray
- For HP/CA burs (ø2.35)

RRP	BUY 1
\$1,955	\$1,662
BUY 2	BUY 3
\$1,564	\$1,466



MODEL **S-Max M65**
ORDER CODE H1008

- Non-Optic
- 1:1 Direct Drive
- Single Spray
- For HP/CA burs (ø2.35)

RRP	BUY 1
\$818	\$695
BUY 2	BUY 3
\$654	\$614



MODEL **FX65m**
ORDER CODE H1014001

- 1:1 Direct Drive
- Clean Head System
- For HP burs (ø2.35)
- Stainless Steel body
- Max speed: 40,000 min⁻¹

RRP	BUY 1
\$461	\$392
BUY 2	BUY 3
\$369	\$346

HANDPIECES

Surgical Handpieces

Implant Handpieces

BUY 1
GET
15% OFF

BUY 2
GET
20% OFF

BUY 3
GET
25% OFF



MODEL **Ti-Max X-SG20L**
ORDER CODE **C1003**

- Optic
- Reducing 20:1
- Titanium
- CA Burs

RRP **\$2,720**
BUY 1 \$2,312
BUY 2 \$2,176 **BUY 3 \$2,040**



MODEL **S-Max SG20**
ORDER CODE **C1010**

- Non Optic
- Reducing 20:1
- Stainless Steel
- CA Burs

RRP **\$1,657**
BUY 1 \$1,408
BUY 2 \$1,326 **BUY 3 \$1,243**



MODEL **Ti-Max Z-SG45L**
ORDER CODE **C1107**

- Optic
- Increasing 1:3
- Triple Spray
- 25mm FG Burs

RRP **\$2,224**
BUY 1 \$1,888
BUY 2 \$1,777 **BUY 3 \$1,666**



MODEL **Ti-Max X-SG93L**
ORDER CODE **C1004**

- Optic
- Increasing 1:3
- Triple Spray
- FG Burs

RRP **\$2,613**
BUY 1 \$2,221
BUY 2 \$2,090 **BUY 3 \$1,960**



MODEL **Ti-Max X-SG65L**
ORDER CODE **H1009**

- Optic
- Direct Drive 1:1
- HP Burs

RRP **\$2,374**
BUY 1 \$2,015
BUY 2 \$1,897 **BUY 3 \$1,778**



MODEL **SGA-E2S**
ORDER CODE **H265**

- Non Optic
- Increasing 1:2
- 65mm HP Burs

RRP **\$2,472**
BUY 1 \$2,101
BUY 2 \$1,978 **BUY 3 \$1,854**



MODEL **SGS-ES**
ORDER CODE **H264**

- Non Optic
- Direct Drive 1:1
- 65mm HP Burs

RRP **\$1,502**
BUY 1 \$1,292
BUY 2 \$1,216 **BUY 3 \$1,140**



MODEL **SGS-E2S**
ORDER CODE **H266**

- Non Optic
- Increasing 1:2
- 20° angulation
- 65mm HP Burs

RRP **\$1,782**
BUY 1 \$1,515
BUY 2 \$1,426 **BUY 3 \$1,337**



MODEL **SGA-ES**
ORDER CODE **H263**

- Non Optic
- Direct Drive 1:1
- 20° angulation
- 65mm HP Burs

RRP **\$2,107**
BUY 1 \$1,791
BUY 2 \$1,686 **BUY 3 \$1,580**


Oral Hygiene Scaler Handpieces



MODEL **VA2-LUX-HP**
ORDER CODE **E351050**

- Optic LED

RRP **\$1,822**
BUY 1 \$1,549 **BUY 2 \$1,457** **BUY 3 \$1,366**



MODEL **VA2-HP**
ORDER CODE **E350050**

- Non-Optic

RRP **\$1,060**
BUY 1 \$901 **BUY 2 \$848** **BUY 3 \$795**

Wide choice of Ultrasonic Scaler Tips

Buy 3 tips, get 1 FREE

MODEL	RRP
Ultrasonic Scaler Tips	From \$128



Scan the QR code for full tip guide

MODEL	ORDER CODE	SALE PRICE
Varios Torque Wrench	Z221076	\$112 RRP-\$132

15% OFF



Upgrade to Optic
or water adjustment
coupling and
get 20% off

20% OFF

OPTIC	MODEL	ORDER CODE	RRP	SALE PRICE	SAVE
•	PTL-CL-LED	P1001600	\$468	\$374	\$94
•	Stainless Steel Body				
•	Titanium coating				
•	Anti-Retraction Valve				
•	PTL-CL-LEDIII	P1001601	\$709	\$567	\$142
•	Water volume adjuster				
•	Stainless Steel Body				
•	Titanium coating				
•	Anti-Retraction Valve				



Varios Combi Pro Basic Kit

2 Maintenance Systems in 1

Varios Combi Pro combines ultrasonic scaling and powder therapy for a more effective and comfortable oral hygiene and biofilm control.

Basic Kit

Model **Varios Combi Pro Basic Kit**
ORDER Code **Y1004194**

- Advanced optic iPiezo scaling for General, Perio and Endo treatment.
- Powder Therapy to clean Supragingival areas efficiently and quickly.
- User friendly handpiece and nozzle for easy access.



MODEL	ORDER CODE	SALE PRICE	SAVE
Varios Combi Pro Basic Kit	Y1004194	\$5,892 RRP \$7,856	\$1,964



Scan the QR code
for more information

Varios 970

Portable/Bench Top

Optimum clinically versatile self-contained portable bench top unit, able to operate using two 400ml solution bottles with independent pumps and an optional water connection.

- iPiezo engine ensures low vibration and stabilised power
- Comfortable, lightweight handpiece
- Easy switch from Endo to Perio to Scaling and General
- LED illumination
- Memory function and LCD display



MODEL	ORDER CODE	SALE PRICE	SAVE
Optic Complete Set	Y1001168	\$3,395 RRP \$4,526	\$1,131
Non-Optic Complete Set	Y1001175	\$2,807 RRP \$3,743	\$936

Prophy-Mate neo

Tooth Polishing System

Prophy-Mate neo, the best polishing tool for beautiful white teeth.

- Removal of plaque and discolouration
- Etching of enamel, dentine and porcelain prior to bonding
- Quick connect to major brands of turbine couplings
- Preparation for pit and fissure sealant
- Sturdy and reliable
- Simple to operate
- Easy maintenance
- Available with 60° or 80° nozzles




MODEL	ORDER CODE	SALE PRICE	SAVE
PMNG-PTL-P for NSK Coupling	Y135029	\$1,272 RRP \$1,696	\$424

Perio-Mate

Biofilm Eraser

Perio-Mate, the best biofilm eraser for gentle periodontal treatment.

- Treatment of periodontal pockets and peri-implantitis
- Eradication of biofilm in periodontal pockets
- Removal of plaque and discolouration
- NSK's innovative nozzle tip
- Precise, flawless powder flow
- Slim nozzle design
- Smooth 360° rotation
- Powder Control Ring



MODEL	ORDER CODE	SALE PRICE	SAVE
PER-PTL-P for NSK Coupling	Y1002654	\$1,272 RRP \$1,696	\$424

Ultimate XL Bench Top Control System

- Control unit • Micromotor Handpiece • Foot Control
- Handpiece Stand • Speed of 1,000-50,000 min⁻¹
- Digital Speed Display • Maintenance Mode Setting
- Auto Cruise Function • Error Code Display Function
- Brushless

MODEL	ORDER CODE	SALE PRICE	SAVE
Ultimate XL Torque micromotor: UMXL-GT (230V)	Y141469	\$2,199 <small>RRP \$2,932</small>	\$733
Ultimate XL Compact micromotor: UMXL-GC (230V)	Y141470	\$2,199 <small>RRP \$2,932</small>	\$733



Ultimate XL-D Bench Top Tower Control System

- Control unit • Micromotor Handpiece • Foot Control
- Handpiece Stand • Speed of 1,000-50,000 min⁻¹
- Digital Speed Display • Maintenance Mode Setting
- Auto Cruise Function • Error Code Display Function
- Available to connect with an air turbine adaptor • Brushless

MODEL	ORDER CODE	SALE PRICE	SAVE
Ultimate XL-D Torque micromotor: UMXL-DT (230V)	Y141504	\$2,199 <small>RRP \$2,932</small>	\$733
Ultimate XL-D Compact micromotor: UMXL-DC (230V)	Y141505	\$2,199 <small>RRP \$2,932</small>	\$733



25% OFF
ULTIMATE XL
RANGE &
VOLVERE i7

Ultimate XL-K Knee Control Top Control System

- Control unit • Micromotor Handpiece • Foot Control
- Handpiece Stand • Speed of 1,000-50,000 min⁻¹
- Digital Speed Display • Maintenance Mode Setting
- Auto Cruise Function • Error Code Display Function
- Available to connect with an air turbine adaptor • Brushless

MODEL	ORDER CODE	SALE PRICE	SAVE
Ultimate XL-K Torque micromotor: UMXL-KT (230V)	Y141495	\$2,199 <small>RRP \$2,932</small>	\$733
Ultimate XL-K Compact micromotor: UMXL-KC (230V)	Y141496	\$2,199 <small>RRP \$2,932</small>	\$733



Ultimate XL-F Floor Type Foot Control System

- Control unit • Micromotor Handpiece • Foot Control
- Handpiece Stand • Speed of 1,000-50,000 min⁻¹
- Auto Cruise Function • Available to connect with an air turbine adaptor • Brushless

MODEL	ORDER CODE	SALE PRICE	SAVE
Ultimate XL-F Torque micromotor: UMXL-FT (230V)	Y141513	\$2,199 <small>RRP \$2,932</small>	\$733
Ultimate XL-F Compact micromotor: UMXL-FC (230V)	Y141514	\$2,199 <small>RRP \$2,932</small>	\$733



VOLVERE i7 Coreless Micromotor

- Light and compact to fit into the smallest space
- Speed 1,000-35,000 min⁻¹ • Microcomputer Control
- Auto Cruise Function • Contents: Control unit, Micromotor Handpiece, Foot control, Handpiece Stand

MODEL	ORDER CODE	SALE PRICE	SAVE
VOLVERE i7 RM	Y1003493	\$1,586 <small>RRP \$2,115</small>	\$529

E-type option available.



MAINTENANCE

iCare C2/C3

Handpiece oiling and maintenance system

- Pre-set solution volume based on handpiece type
- Thorough cleaning and lubrication with simple operation
- Chuck cleaning and lubricating function
- Purging mode removes excess oil
- Integrated 1.2 litre oil tank reduces burden of refill

25% OFF
WITH Q2 PROMO



MODEL	ORDER CODE	SALE PRICE	SAVE
iCare C2	Y1002787	\$4,117 RRP-\$5,489	\$1,372
<ul style="list-style-type: none"> • Standard System Complete • 2 rotation ports • 2 non-rotation ports • 230V • Requires 2 adaptors 			

iCare C3	Y1002788	\$4,397 RRP-\$5,863	\$1,466
<ul style="list-style-type: none"> • Standard System Complete • 3 rotation ports • 1 non-rotation port • 230V • Requires 3 adaptors 			

MODEL	ORDER CODE	SALE PRICE	SAVE
Adaptor for NSK High Speed Handpieces that fit NSK FlexiQuik type coupling	Z257010	\$179 RRP-\$238	\$60
Adaptor for W&H High Speed to fit W&H Roto Quick Coupling	Z257021	\$179 RRP-\$238	\$60
Adaptor for Bien Air High and Low Speed that fit Bien Air Unifix Coupling	Z257022	\$179 RRP-\$238	\$60
Adaptor for Sirona High and Low Speed that fit Sirona Quick Coupling	Z257023	\$179 RRP-\$238	\$60
Adaptor For Endo-Mate DT/TC, Taskal Wizards Heads	Z257042	\$179 RRP-\$238	\$60
Adaptor to fit Morita High Speed Handpiece	Z257048	\$179 RRP-\$238	\$60

Pana Spray offer

Buy 1 box of 6,
get 1 box of 6
Pana Spray Plus FREE

MODEL **Pana Spray Plus Pack**
ORDER CODE Y900630

- For High and Low speed handpieces, and air motors
- 500 mL • Alcohol-contained lubricant
- Pack of 6 cans

SALE PRICE
\$269
SAVE
\$269
RRP-\$538



**BUY 1 BOX
GET 1 BOX
FREE**

ENDODONTICS

Simple and precise root canal treatment

Endodontic Motors

ENDO-MATE TC2

Model **ENDO-MATE TC2 (230V)**
MP-F16R head

ORDER Code Y1001027

- Large LCD display with 5 programs for different file systems
- Supports most major brands of Ni-Ti files
- Torque control facility ranging from 0.1 to 4.0 NCM depending on head:gear ratio



SALE PRICE
\$1,697
SAVE
\$566
RRP-\$2,263

ENDO-MATE AT

Model **ENDO-MATE AT (230V)**
ORDER Code Y1002190

- Fits all major brands of Files
- Torque control, with integrated auto-reverse function
- Compact stylish control/monitor pad for simple and user friendly operation



SALE PRICE
\$1,384
SAVE
\$461
RRP-\$1,845

Apex Locator

iPex II

Model **iPex II Complete set**
ORDER Code Y1002208

- Accurate measurement of any root canal
- SmartLogic technology
- Operates in all conditions, including dry, wet and bloody canals



SALE PRICE
\$1,272
SAVE
\$424
RRP-\$1,696

25% OFF

"The NSK Endomate TC2 has been my motor of choice for the past 10 years. It is a well balanced, cordless motor that delivers consistent torque and speed for the endodontic treatment."



Dr Morgan Wu
Specialist Endodontist
Bds (Otago), Dcd (Endo)(Melb)

NSK Oceania Pty Ltd www.nskoceania.com.au

Product Information | Technical Support

Freephone: 1300 44 33 21 (Australia) Email: info@nskoceania.com.au

All offers as listed. Valid until 30 June 2022. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.