

Special Offers

January - March 2023

Z95L
Contra-angle



Z900L
Air Turbine

DualDefense™

Proprietary Infection Control by NSK

NSK instruments are designed with DualDefense infection protection.

The Clean Head System prevents suck-back into the head and the exhaust line


and the Anti-Retraction valve prevents contamination of the water supply.

Clean Head System

Anti-suck back system prevent cross contamination of the air lines

Anti-Retraction Valve

Anti-Retraction systems prevents cross contamination of the water lines in High-Speed Air Turbines and Electric Micromotors.



All NSK handpieces are equipped with a **Clean Head System** that prevents suck-back.

All coupling and micromotors have integrated **Anti-Retraction systems**.



The efficacy of the DualDefense system has been tested and validated by well known universities, such as the University of Washington, Umeå University in Sweden and National University of Singapore.

Handpieces

Ti-Max Z

Titanium Turbine



Power 26W

- Clean Head System (Anti-Retraction Valve)
- Torque: 26W (model Z900L)
- Smallest head : $\phi 12.5 \times H 13.1$ mm
- Speed: 320,000~400,000 min⁻¹
- Titanium body with DURAGRIP
- Quattro Water Spray
- Push Button chuck
- Ceramic bearings
- Microfilter
- **2 Years warranty**

NSK connection

MODEL **Z900L**
ORDER CODE **P1111**
RRP **\$1,981**
SALE PRICE **\$1,525**

KaVo® Multiflex® connection

MODEL **Z900KL**
ORDER CODE **P1113**
RRP **\$2,074**
SALE PRICE **\$1,659**



Ti-Max Z

Titanium Contra Angle



- Clean Head System (Anti-Retraction)
- Titanium body with scratch-resistant DURAGRIP
- Cellular Glass Optics
- Push Button Chuck
- Microfilter
- **2 Years warranty**

Switch Function

- Ti-Max Z95L and Ti-Max Z45L only
- Switch between Mist Spray and Jet Spray
- Mist Spray (Water and air)
- Jet Spray (Water only)
- No equipment modification required

MODEL **Z24L**
ORDER CODE **C1126001**

- Optic
- 1:1 Direct Drive
- Mini Head
- Single spray
- For CA burs ($\phi 2.35$)

RRP **\$1,491**
SALE PRICE **\$1,193**



MODEL **Z25L**
ORDER CODE **C1038**

- Optic
- 1:1 Direct Drive
- Single spray
- For CA burs ($\phi 2.35$)

RRP **\$1,491**
SALE PRICE **\$1,193**



MODEL **Z95L**
ORDER CODE **C1034SW**

- Optic
- 1:5 Speed Increasing
- Switch Spray function
- For FG Burs
- **3-Year Warranty**

RRP **\$2,121**
SALE PRICE **\$1,697**



MODEL **Z45L**
ORDER CODE **C1064**

- Optic
- 1:4.2 Speed Increasing
- Switch Spray function
- 45° Head
- FG Burs

RRP **\$2,263**
SALE PRICE **\$1,810**



Couplings



Optic NSK connection
MODEL **PTL-CL-LED**
ORDER CODE **P1001600**

- Anti-retraction value

RRP **\$468**
SALE PRICE **\$374**



Optic NSK connection
MODEL **PTL-CL-LEDIII**
ORDER CODE **P1001601**

- Anti-retraction value
- Water volume adjuster

RRP **\$709**
SALE PRICE **\$567**



Optic KaVo connection
MODEL **KCL-LED**
ORDER CODE **P1005600**

- Anti-retraction value
- Water volume adjuster

RRP **\$795**
SALE PRICE **\$636**

Electric Micromotor

NLZ E

Electric Micromotor with Endodontic Function

Simple installation

- All clinical application use
- Reciprocating and rotational Endo Function
- Contra-Check function
- High speed rotation and torque
- Flexible speed variation
- No maintenance
- ISO Type E easy to connect with all handpieces and contra-angles (NSK, KaVo®, Sirona®, W&H®, Bien-Air®)
- Simple installation to existing high speed line

MODEL **NLZ E**
ORDER CODE **Y1500361**
RRP **\$5,942**
SALE PRICE **\$4,754**



Surgical

Osseo 100

Osseointegration measuring device

Removes doubt

- Osseo 100 measures implant stability and osseointegration to enhance decisions about when to load the implant. Especially important when using protocols with shorter treatment time and treating higher risk patients.
- By mounting a MultiTeg™ the measurement is made in a second. A numerical result appears on an LED screen and that indicates how well the implant has integrated.

MultiTegs sold separately – for full list scan the QR code.



MODEL **Osseo 100**
ORDER CODE **Y1004175**
RRP **\$3,806**
SALE PRICE **\$3,045**



iSD900

Cordless Torque Screwdriver

Faster and safer

- iSD900 cordless screwdriver with Torque Calibration System (TCS) for prosthodontic screw fixing procedures. iSD900 will guarantee precise torque control due to the unique NSK torque calibration system.
- Compatible with all major implant systems
- With torque calibration (TCS)
- Precise torque between 10 and 40 Ncm
- Speeds of 15, 20 to 25 min⁻¹

MODEL **iSD900**
ORDER CODE **Y1001358**
RRP **\$3,604**
SALE PRICE **\$2,883**



Varios 970 iPiezo engine® Multifunctional Ultrasonic System

Portable

- All-in-one unit for Periodontic, Endodontic, Scaling, MI applications and Implant maintenance.
- Convenient large capacity bottles of 400ml of the same solution each or different solutions with independent pumps.
- The iPiezo engine assures stable output by automatically adjusting to the optimum frequency depending on the loading of the tip.
- Includes tips G4, G6 and G8
- Scan the QR code to view full tip range.



MODEL **Varios 970 Optic**
ORDER CODE **Y1001168**
RRP **\$4,526**
SALE PRICE **\$3,621**



Varios Combi Pro

Two systems in one

Combining ultrasonic scaling and powder therapy for more effective and comfortable oral hygiene and biofilm control.

Bundle offer

Includes;

- Varios Combi Pro Unit
- PROPHY-Chamber for Supragingival air polishing
- x 3 PROPHY-Handpieces with 60° Nozzle for supra gingival cleaning
- x4 300g bottles of FLASH pearl supragingival cleaning powder
- PERIO-Chamber for subgingival air polishing
- x2 PERIO-Handpieces
- x 2 boxes (40 pcs) Sterile disposable nozzle tips
- x2 160g bottles of Perio-Mate powder
- x2 VA2LUX-Optic Scaler Handpieces
- 3 Scaling tips (G4, G6, G8)and sterilisation case

Plus Free iCart V

ORDER CODE **VCPBUNDLEPKG**
WORTH **\$12,751**
SAVE **\$2,536**



"Stains have never been so easy to remove! I love how I can customize the power and water easily during treatment and between patients, allowing me to provide my patients with a pleasant cleaning appointment, best suited for their needs. Highly recommend the Varios Combi Pro, can't go without it now!"



Tessa Cowburn
Oral Health Therapist



NSK Products are precision manufactured in Japan with 100% of NSK parts manufactured in-house. It is with this quality in mind that NSK products come with guarantees of the highest standard. Registering your NSK products today for Warranty will ensure total peace of mind. With NSK, you are supported for the life of your NSK products.



To register your product for warranty visit **NSKtech.com.au** or simply scan the QR Code

**Your local NSK
Business Development
and Education
Specialists**

WA
Stuart Blair
0437 407 668
stuart@nskoceania.com.au

Queensland
Sam Vaizie
0439 260 006
sam@nskoceania.com.au

NSW North
Theresa Sparkes
0419 251 216
theresa@nskoceania.com.au

VIC, SA & NT
Mariam Dalzotto
0427 833 684
mariam@nskoceania.com.au

NSW South and ACT
Janette Carroll
0439 199 979
janette@nskoceania.com.au

Tasmania
Stuart Blair
0437 407 668
stuart@nskoceania.com.au

CONNECT WITH US



NSK Oceania Pty Ltd www.australia.nsk-dental.com

Product Information | Technical Support

Freephone: 1300 44 33 21 (Australia) Email: info@nskoceania.com.au

Offers available to listed items only. Valid until 31 March 2023. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.