

BLACK FRIDAY PROMO OFFERS

UPTO
25%
OFF RRP

Ti-Max Handpieces



MODEL Ti-Max Z990L
ORDER CODE PA23730001

- Standard Head
- 44W*
- 2-year warranty

RRP \$1981

Sale Price
\$1,486



MODEL Ti-Max Z890L
ORDER CODE PA23720001

- Miniature Head
- 41W*
- 2-year warranty

RRP \$1981

Sale Price
\$1,486



MODEL Ti-Max ZMicro
ORDER CODE PA2357

- Micro Head
- 20W
- 2-year warranty

RRP \$1381

Sale Price
\$1,036



MODEL Ti-Max X450L
ORDER CODE P1078

- Rear Exhausting Head
- 21W
- 2-year warranty

RRP \$1677

Sale Price
\$1,258



MODEL Ti-Max Z95L
ORDER CODE C10345W

- Mist or Jet Spray Switch
- 1:5 Increasing
- 3-year Warranty

RRP \$2121

Sale Price
\$1,591



MODEL Ti-Max Z25L
ORDER CODE C1038

- Single Water Spray
- 1:1 Direct Drive
- 2-year Warranty

RRP \$1491

Sale Price
\$1,118



MODEL Ti-Max Z24L
ORDER CODE C1126001

- Miniature Head
- Single Water Spray
- 2-year warranty

RRP \$1491

Sale Price
\$1,118



MODEL Ti-Max Z45L
ORDER CODE C1064

- Mist or Jet Spray Switch
- 1:4.2 Increasing
- 2-year warranty

RRP \$2263

Sale Price
\$1,697

S-Max Handpieces




MODEL S-Max M900L
ORDER CODE P1254

- Standard Head
- 26W
- 1-year warranty

RRP \$1381

Sale Price
\$1,036



MODEL S-Max M800L
ORDER CODE P1253

- Miniature Head
- 23W
- 1-year warranty

RRP \$1381

Sale Price
\$1,036



MODEL S-Max Pico
ORDER CODE P1140

- Micro Head
- 9W
- 1-year warranty

RRP \$1381

Sale Price
\$1,036



MODEL S-Max M95L
ORDER CODE C1023

- Quattro Spray
- 1:5 Increasing
- 1-year Warranty

RRP \$1763

Sale Price
\$1,322



MODEL S-Max M25L
ORDER CODE C1024

- Single Water Spray
- 1:1 Direct Drive
- 1-year Warranty

RRP \$1353

Sale Price
\$1,015



MODEL S-Max M25
ORDER CODE C1027

- Non-Optic
- Single Water Spray
- 1-year warranty

RRP \$949

Sale Price
\$712



MODEL S-Max M65
ORDER CODE H1008

- HP Burs
- Single Water Spray
- 1-year warranty

RRP \$818

Sale Price
\$614

High-Speed
Connections
available for KaVo,
W&H, Sirona,
Bien-Air

iCare

NSK's iCare simplifies the maintenance of dental handpieces. Automatically lubricate up to four handpieces simultaneously, ensuring optimal performance with minimal effort. Compact and stylish, the iCare fits seamlessly into any dental practice, with easy setup and maintenance. Available in 2 configurations and adaptable to various handpieces, the iCare is a versatile and essential unit to increase the lifespan of your handpieces.

Specifications

- Power Supply AC 100-240 V 50/60Hz • Air Pressure 0.5-0.6 MPa (5.0-6.0 kgf/cm²) • Weight 6.5 kg (C2 Type / C3 Type) • Dimension W 280 x D 240 x H 350 mm • Solution Tank Capacity 1.2 L

Contents

- iCare Unit • Power Cord • Chuck Lubrication Nozzle • Maintenance Oil for iCare (1 Litre) • Mist Filter Set • Oil Absorber Sheet • Air Tube
- * Connects only to Midwest 4-hole type
- * Air piping is necessary (5.0-6.0 kgf/cm²)



MODEL **iCare C2**

ORDER CODE Y1002787

RRP \$5489

SALE PRICE

\$4,117

MODEL **iCare C3**

ORDER CODE Y1002788

RRP \$5863

SALE PRICE

\$4,397

PANA SPRAY Plus

PANA SPRAY Plus is a specially formulated, highly effective general handpiece and air motor lubricant. Frequent application of PANA SPRAY Plus will assure handpiece and air motor performance and durability. Apply PANA SPRAY Plus every time before autoclaving any handpiece or air motor.

- For High & Low speed handpieces, and air motors Alcohol-contained lubricant Contains 480 mL Each

MODEL **PANA SPRAY Plus (Box of 6)**

ORDER CODE Y900630G

RRP \$269

SALE PRICE

\$242



iClave

The iClave Plus by NSK offers advanced sterilization with an 18-liter capacity, 50% larger than conventional autoclaves. Its user-friendly panel features color-coded indicators for reliable operation. The Adaptive Heat System ensures even heating with minimal loss. All sterilization cycles are recorded on a USB drive for easy access. The vacuum sterilization system, combined with a bacterial filter, prevents re-contamination. With multiple programs for different needs, the iClave Plus ensures thorough and efficient sterilization for all dental instruments.

Specifications

- External Dimensions W 445 x D 532 x H 428 mm
- Chamber Dimensions Ø240 x 384 mm
- Chamber Capacity 18 Litre
- Net Weight 55 kg
- Maximum Power Consumption 1,900 W
- Supply Voltage CE 230 V 50Hz
- Air Expulsion System Vacuum Pump 1, 3, 4
- Vacuum Max Load Solid 4 kg
- Max Load Porous 1.5 kg

MODEL **iClave Plus**

ORDER CODE Y1003077

RRP \$10,692

SALE PRICE

\$9,623



NSK Create it

Endodontics



Endodontic Handpieces



MODEL **NLZ Endo**
ORDER CODE **C1130**

- Non-Optic
- 6:1 Reduction
- Rotational files

RRP \$893

Sale Price
\$804



MODEL **FX75m**
ORDER CODE **C1055001**

- Non-Optic
- 16:1 Reduction
- Rotational files

RRP \$788

Sale Price
\$709



MODEL **MP-F16R**
ORDER CODE **Y110096**

- Fits ENDO-MATE Systems
- 16:1 Reduction
- Rotational files

RRP \$763

Sale Price
\$687



MODEL **MP-F20R**
ORDER CODE **Y110044**

- Fits ENDO-MATE Systems
- 20:1 Reduction
- Rotational files

RRP \$763

Sale Price
\$687



MODEL **MPAS-F20R**
ORDER CODE **Y110046**

- Fits ENDO-MATE Systems
- Apex Locator Connection
- 20:1 Reduction
- Rotational files

RRP \$1043

Sale Price
\$939



MODEL **MPAS-F16R**
ORDER CODE **Y110098**

- Fits ENDO-MATE Systems
- Apex Locator Connection
- 16:1 Reduction
- Rotational files

RRP \$1043

Sale Price
\$939



MODEL **MPA-F20R**
ORDER CODE **Y110045**

- Fits ENDO-MATE Systems
- Apex Locator Connection
- 20:1 Reduction
- Rotational files

RRP \$1043

Sale Price
\$939



MODEL **MPA-F16R**
ORDER CODE **Y110097**

- Fits ENDO-MATE Systems
- Apex Locator Connection
- 16:1 Reduction
- Rotational files

RRP \$1043

Sale Price
\$939

iPexII

The Advanced Apex Locator features NSK's SmartLogic technology, providing precise root apex detection across various tooth forms. It accurately displays the file tip's location on its LCD panel, based on extensive clinical trials, ensuring reliable endodontic measurements.

- Control Unit • iPexII Probe (1.8 m) • iPexII File Clip (3 pcs.) • Lip Hook (3 pcs.) • AAA Manganese Dry Batteries (3 pcs.) • iPexII Tester

MODEL **iPexII**
ORDER CODE **Y1002208**
RRP \$1696

SALE PRICE
\$1,272



ENDO-MATE DT2

The ENDO-MATE DT2, compatible with major Ni-Ti file brands, features a pre-programming function for efficient use of various file protocols. A lightweight, compact control unit with a rechargeable battery for enhanced mobility. The slim, lightweight handpiece ensures comfortable handling and high performance.

- Control unit • Motor and cord • MP-F20R head (20:1)
- Handpiece stand • AC adaptor

Specifications

- Power Supply 230V 50/60 Hz and Battery pack
- Number of Programs 9
- Max Torque 7 Ncm (when using 20:1 head)
- Speed range 100-13,000 min-1 (With 20:1, 1:1 geared heads)
- Charge/Continuous use 5 h / 2 h (depending on workload)
- Dimensions W 93 x D 147 x H 100 (mm)
- Weight Unit 440 g / Motor handpiece and Cord 95 g

MODEL **Endo-Mate DT2**

ORDER CODE Y1004230

RRP \$2348

SALE PRICE **\$1,761**



ENDO-MATE AT

The ENDO-MATE AT endodontic motor is compatible with all major Ni-Ti file brands and offers adjustable torque and speed settings. Its compact control pad features a large, easy-to-use screen and an auto-reverse function. The lightweight handpiece enhances handling and visibility, and it includes a convenient on/off switch for streamlined operation.

- Control unit • Motor and cord • MP-F20R head (20:1)
- Handpiece stand • AC adaptor

Specifications

- Power Supply AC 230 V 50/60 Hz
- Number of Programs 9
- Max Torque 7 Ncm (When Using 20:1 Head)
- Speed Range 100-13,000 min-1
- Dimensions W 268 x D 230 x H 103 mm
- Weight 342 g (Including Motor handpiece & cord)

MODEL **Endo-Mate AT**

ORDER CODE Y1002190

RRP \$1845

SALE PRICE **\$1,384**



ENDO-MATE TC2

The Cordless ENDO-MATE TC2 Endodontic Handpiece enhances safety and efficiency with advanced torque control, auto-reverse, and a slim design for improved visibility. It features programmable settings, dual alarm sounds, and compatibility with most Ni-Ti files. An integrated LCD panel and the ability to combine with an Apex locator facilitate precise treatments. The device includes a quick 90-minute non-contact charger for fast recharging.

Contents

- ENDO-MATE TC2: • MP-F16R Head • ENDO-MATE TC2 Motor • Battery Charger
- ENDO-MATE TC2 MPA : • MPA-F16R Head • ENDO-MATE TC2 Motor • Battery Charger

Specifications

- Power Supply 230V 50/60 Hz and Battery pack
- Number of Programs 9
- Max Torque 7 Ncm (when using 20:1 head)
- Speed range 100-13,000 min-1 (With 20:1, 1:1 geared heads)
- Charge/Continuous use 5 h / 2 h (depending on workload)
- Dimensions W 93 x D 147 x H 100 (mm)
- Weight Unit 440 g / Motor handpiece and Cord 95 g

MODEL **Endo-Mate TC2**

ORDER CODE Y1001027

RRP \$2263

SALE PRICE **\$1,697**

MODEL **Endo-Mate TC2 with MPA**

ORDER CODE Y1001028

RRP \$2716

SALE PRICE **\$2,037**



NLZ-E

NSK's "NLZ-E" system offers a seamless transition from air-driven to electric handpieces, enhancing safety and flexibility in dental practices. Featuring an advanced micromotor technology, the NLZ-E broadens the scope of dental treatments with its innovative design tailored for rotational and reciprocating endodontic applications. This system ensures improved comfort and safety, requiring only a minimal update to existing office equipment to significantly upgrade the treatment environment.

Contents

- NLZ E U (Control Unit & Main Unit) • NE180 (AC Adaptor)
- NLZ (Motor) • NLZ CD (Motor Cord) • NLZ U Bracket

Specifications

- Dimension Control Unit W 76 x D 87.1 x H 63.4 mm
- Dimension Main Unit W 146 x D 75 x H 48 mm
- Power Supply 230 V 50/60 Hz*
- Programs General Application Mode / 6 Programs Rotary Endo Mode / 6 Programs
- Endodontics Torque 0.3-6.0 Ncm (Depending on Gear Ratio)
- Motor Speed 100-40,000 min-1

MODEL **NLZ-E**

ORDER CODE Y1500361

RRP \$5942

SALE PRICE **\$4,456**



NSK Create it



Surgical Handpieces



MODEL **Ti-Max X-SG20L**
ORDER CODE **C1003**

- Optic
- 20:1 Reduction
- CA Burs

RRP \$2,720

Sale Price
\$2,040



MODEL **S-Max SG20**
ORDER CODE **C1010**

- Non-Optic
- 20:1 Reduction
- CA Burs

RRP \$1,657

Sale Price
\$1,243



MODEL **Ti-Max Z-SG45L**
ORDER CODE **C1107**

- Optic
- 1:3 Increasing
- 25mm FG Burs

RRP \$2,221

Sale Price
\$1,666



MODEL **Ti-Max X-SG65L**
ORDER CODE **H1009**

- Optic
- 1:1 Direct Drive
- HP Burs

RRP \$2,371

Sale Price
\$1,778



MODEL **SGA-ES**
ORDER CODE **H263**

- Non-Optic
- 1:1 Direct Drive
- 65mm HP Burs

RRP \$2,107

Sale Price
\$1,580



MODEL **SGA-E2S**
ORDER CODE **H265**

- Non-Optic
- 1:2 Increasing
- 65mm HP Burs

RRP \$2,472

Sale Price
\$1,854



MODEL **SGS-ES**
ORDER CODE **H264**

- Non-Optic
- 1:1 Direct Drive
- 65mm HP Burs

RRP \$1,520

Sale Price
\$1,140



MODEL **SGS-E2S**
ORDER CODE **H266**

- Non-Optic
- 1:2 Increasing
- 65mm HP Burs

RRP \$1,782

Sale Price
\$1,337



MODEL **SGX-E20R**
ORDER CODE **HA1200**

- Non-Optic
- 20:1 Reduction
- CA Burs

RRP \$3,388

Sale Price
\$2,541



iSD900

The iSD900 cordless prosthodontic screwdriver features a Torque Calibration System (TCS) for precise control during screw fixing procedures. Designed for easy use, ensuring precise and risk-free screw tightening, which enhances both safety and productivity in clinical settings.

Contents

- iSD900 Motor • iSD-HP • Quick Charger for iSD900
- Torque Calibrator • On/Off Switch Lever

MODEL **iSD900**
ORDER CODE **Y1001358**
RRP \$3,604

Sale Price
\$2,703



Osseo 100+

The Osseo 100 device measures implant stability and osseointegration, facilitating timely decisions on implant loading. Ideal for shorter treatment schedules or high-risk patients. The non-invasive ISQ measurement process aids in planning the implant loading period to minimize failure risks. Links with Surgic Pro2

Contents

- Osseo 100+ Instrument • MultiTipeg Driver
- Mains adapter and plugs

MODEL **Osseo100+**
ORDER CODE **Y1004176**
RRP \$3,806

Sale Price
\$2,854

MODEL **Osseo100**
ORDER CODE **Y1004175**
RRP \$3,267

Sale Price
\$2,450

MultiTipeg sold separately

VarioSurg3

VarioSurg3 Ultrasonic bone surgery system provides safe and precise cutting of bone, without damaging soft tissue. A variety of tips makes procedures easier and stress-free to help provide more precise and efficient treatment.

Contents

- Control Unit • LED Handpiece with 2m cord • FC-78 Foot Control • Sterilisation Cassette • Handpiece Stand • Irrigation tube (5 pcs.) and other accessories • Basic H-S Kit (H-SG1, SG3, SG5 SG6D, SG7D, SG11 and Tip holder)

MODEL **VarioSurg3 complete set W/O foot control**

ORDER CODE Y1002248

RRP \$10,060

SALE PRICE \$7,545

MODEL **VarioSurg3 complete set**

ORDER CODE Y1002726

RRP \$12,468

SALE PRICE \$9,351



Surgic Pro2

The Surgic Pro2 improves operability, safety, and connectivity. It features a wireless foot control, Advanced Handpiece Calibration for precise torque settings and integrates with the Osseo 100+ device for real-time, contactless monitoring of implant stability via Bluetooth.

Contents

- Control Unit • Micromotor • 20:1 Reduction Handpiece • FOOT CONTROL WIRELESS • Irrigation tube (pack of 3) and other accessories

MODEL **Surgic Pro2 OPTIC**

ORDER CODE Y1500511

RRP \$8,726

SALE PRICE \$6,544

MODEL **Surgic Pro NON-OPTIC**

ORDER CODE Y1500512

RRP \$7,479

SALE PRICE \$5,609



VarioSurg3 Tip Kits

MODEL **Sinus Lift Kit**

- Tip Holder • SG1, SG3, SG6D, SG9, SG10, SG11

ORDER CODE Y900689

RRP \$2,096

SALE PRICE \$1,886


MODEL **Basic H-S Kit**

- Tip Holder • H-SG1, SG3, SG5, SG6D, SG7D, SG11

ORDER CODE Y1002775

RRP \$2,183

SALE PRICE \$1,965



iCart Duo

The NSK iCart Duo holds both the VarioSurg3 & Surgic Pro2. Use one foot control and link together with the Link Module.

MODEL **iCart Duo**

ORDER CODE S9090

RRP \$2,354

SALE PRICE

\$1,766

MODEL **Link Module**

ORDER CODE Z1402001

RRP \$331

SALE PRICE

\$298



Increase your workflow with 10% off additional Surgic Pro & VarioSurg3 cables

| MODEL | ORDER CODE | SALE PRICE |
|----------------------------------------------------------------------|------------|-------------------------------|
| VS3-LED-HPSC Optic scaler with 2m cord for VarioSurg3 | E1133 | \$4,965 RRP \$5,517 |
| SGL70M Optic micromotor with 2m cord for Surgic Pro | E1023001 | \$4,018 RRP \$4,464 |
| SG70M Non-Optic micromotor with 2m cord for Surgic Pro | E1025001 | \$3,193 RRP \$3,548 |
| SGL80M Optic micromotor with 2m cord for Surgic Pro2 | Y1004211 | \$4,018 RRP \$4,464 |
| SG80M Non-Optic micromotor with 2m cord for Surgic Pro2 | Y1004212 | \$3,193 RRP \$3,548 |





Prophy Handpieces



MODEL **FX57m**
ORDER CODE **C1054001**

- Non-Optic
- 4:1 Reduction
- Screw in cups

RRP \$508

Sale Price
\$457



MODEL **FX25m**
ORDER CODE **C1052001**

- Non-optic
- 1:1 Direct Drive
- Push Button CA burs

RRP \$462

Sale Price
\$416



MODEL **Ti-Max X57L**
ORDER CODE **C609**

- Optic
- 16:1 Reduction
- Screw in cups

RRP \$1601

Sale Price
\$1,201



MODEL **Ti-Max X55L**
ORDER CODE **C608**

- Optic
- 4:1 Reduction
- IPR Profin/Eva Tips

RRP \$2017

Sale Price
\$1,513



MODEL **EC-30M**
ORDER CODE **C100007**

- Non-Optic
- 1:1 Direct Drive
- Contra-Angle Shank

RRP \$190

Sale Price
\$171



MODEL **ER4M**
ORDER CODE **C634002**

- Non-Optic
- 4:1 Reduction
- Contra-Angle Shank

RRP \$396

Sale Price
\$356



MODEL **ARY (S)**
ORDER CODE **C187**

- Non-optic
- Max Speed 5,000 min
- Screw-in cups

RRP \$113

Sale Price
\$102



MODEL **NAC (Y)**
ORDER CODE **C032**

- Non-optic
- Max Speed 10,000 min
- Latch type CA burs

RRP \$126

Sale Price
\$113

Prophy-Mate Neo

Air Powered Tooth Polishing System

Prophy-Mate Neo combines form and functionality for user convenience. Its excellent weight balance and lightweight, compact powder chamber ensure a comfortable grip. The connection swivels 360°, making it easy to use. This effortless operation reduces the time needed for treating plaque and stains. Available nozzles include 60° and 80°.

- Handpiece with 60° Nozzle Head
- Powder Case
- Replacement Nozzle Head 80°
- 5 pcs. x 15 g Sachets of FLASH pearl Cleaning Powder
- Maintenance Kit



MODEL **PMNG-PTL-P**

ORDER CODE **Y135029**

RRP \$1696

SALE PRICE **\$1,272**

Perio-Mate

Biofilm Eraser

Perio-Mate is used in the treatment of periodontal pockets and peri-implantitis. The Perio-Mate's ultrafine water spray and glycine-based Perio-Mate Powder gently eradicate biofilm in periodontal pockets in just 20 seconds. A safe and supremely efficient spray pressure method that enables fast, precise removal of plaque biofilm from subgingival areas

- Handpiece
- Powder Case
- Nozzle Tip (40 pcs.)
- Maintenance Kit
- Perio-Mate Powder (160 g) x 1



MODEL **PER-PTL-P**

ORDER CODE **Y1002654**

RRP \$1696

SALE PRICE **\$1,272**

Varios Combi Pro

Ultrasonic scaling X Air Polishing

NSK has combined ultrasonic treatment and powder therapy into a single professional treatment unit that streamlines procedures and saves time and space. The detachable chamber and handpiece make it possible to easily switch between supragingival and subgingival applications and make the unit extremely user and maintenance friendly.

Basic Set Contents

- PROPHY-Handpiece with 60° Nozzle
- PROPHY-Hose
- PROPHY-CHAMBER
- FLASH pearl (5 pcs. x 15 g sachets)
- Maintenance kit
- VA2LUX-Handpiece
- VA2LUX-Hose
- 3 Scaling tips
- Sterilization case
- Irrigation bottle (400 mL)
- Foot control

Bundle Package Contents

- Varios Combi Pro Unit
- PROPHY-CHAMBER for supragingival air polishing
- x3 Prophy Handpieces with 60° Nozzle for supragingival cleaning
- x4 300g bottles of FLASH pearl supragingival cleaning powder
- PERIO-CHAMBER for subgingival air polishing
- x2 Perio Handpieces
- x2 boxes (40 pcs) Sterile disposable nozzle tip
- x2 160g bottles of Perio-Mate Powder
- x2 Optic scaler handpieces
- x3 Scaling tips (G4, G6, G8) & Sterilization case

- Power supply : AC 100-240 V 50/60 Hz
- Frequency : 28-32 kHz
- Dimensions : W 225 x D 290 x H 162 (mm) (including chamber)
- Unit weight : 2.6 kg (Except attachment)
- Bottle volume : 400 mL



MODEL **VCP Bundle Package**

ORDER CODE VCPBUNDLEPKG

RRP \$12751

Offer

Free iCartV

MODEL **Varios Combi Pro Basic Set**

ORDER CODE Y1004194

RRP \$7856

SALE PRICE

\$5,892

Varios 970

Portable Ultrasonic Scaler

The functionality required for ultrasonic applications is available in the most-advanced model in our compact series, the Varios 970. Our high-capacity 2 bottle system allows users to focus on safe treatment without worrying about water irrigation. Simple daily maintenance with auto rinse function and seamless body design.

Contents • Control unit • Handpiece • Handpiece cord • 3 Scaling tips (G4, G6, G8) • Foot control (FC-70) • 2 x Irrigation bottles • 3 x Tip wrench with torque limiter • AC power cord • Sterilization case

- Slimmest and lightest* handpiece with LED illumination
- iPiezo engine ensures stable power output and performance
- Twin LED lights through twin optical glass rods
- Optic and non-optic handpieces available
- Extensive tip range for a wide range of applications

MODEL **Varios 970**


ORDER CODE Y1001168

RRP \$4526

SALE PRICE

\$3,394

View our complete tip range



VIVA ace

Mobile Dental Unit

A complete treatment centre in one compact unit. Complete with Suction, Ultrasonic scaler, Multi functional micromotor and 3&1 Syringe, the VIVAace allows more comfortable and advanced treatment, outside of the dental centre.

- E-Type motor, ultrasonic scaler, 3 in 1 syringe and suction combined in a single compact body including storage space for all instruments.
- Lightweight at 8.6 kg
- Dual foot controls to operate micromotor, ultrasonic scaler and suction
- Full endodontic functionality can be performed with (1:1) contra-angle.

MODEL **VIVAace Complete Set**

ORDER CODE Y1004228

RRP \$21060

SALE PRICE

\$17,901

NSK Create it



Laboratory

Straight Handpieces



MODEL **Ti-Max X65L**
ORDER CODE **H273**

- Optic
- Single Water Spray
- HP Burs

RRP \$1955

Sale Price
\$1,466



MODEL **S-Max M65**
ORDER CODE **H1008**

- Non-Optic
- Single Water Spray
- HP Burs

RRP \$818

Sale Price
\$613



MODEL **FX65m**
ORDER CODE **H1014001**

- Non-Optic
- No Water Spray
- HP Burs

RRP \$461

Sale Price
\$345

VOLVERE i7

The VOLVERE i7 by NSK offers dental technicians a compact, stylish, and efficient lab engine. Designed to fit into tight spaces. Features include an auto-cruise function to reduce foot fatigue, microcomputer control for smooth operation, and a powerful, low-vibration micromotor with a dustproof mechanism.

Contents

- Standard Micromotor System • VOLVERE i7 Control Unit • Motor Handpiece • Foot Control • Handpiece Stand
- E-Type Micromotor System • VOLVERE i7 • Control Unit E-Type Micromotor • Foot Control • Handpiece Stand

Specifications

• Control Unit • Power Supply AC 230 V 50/60 Hz • Dimensions W 69 x D 185 x H 167 mm • Weight 930g

• Standard Type • Motor Handpiece Speed 1,000-35,000 min • Max. Torque 4.1 Ncm • Dimensions L 145 x Ø24.5 mm • Weight 192 g (without Cord) • Cord Length 1.2 m

• E-Type Motor Handpiece Speed 1,000-35,000 min • Max. Torque 4.1 Ncm • Dimensions L 109.3 x Ø24.5 mm • Weight 118 g (without Cord) • Cord Length 1.2 m



MODEL **VOLVERE i7 Standard**
ORDER CODE **Y1003493**

RRP \$2115

SALE PRICE
\$1,586



MODEL **VOLVERE i7 E-Type**
ORDER CODE **Y1003494**

RRP \$1741

SALE PRICE
\$1,306

ULTIMATE XL Series

The Ultimate XL laboratory micromotor offers unprecedented smoothness and power with its 180° vector control system, ensuring efficient, stress-free work. Its ergonomic handpiece design reduces hand fatigue, making it ideal for delicate tasks. Features include low noise and vibration, a compact, lightweight body, digital speed display, and an adaptable turbine for versatile control.



NSK Ultimate XL Bench Top

Contents

- Control Unit • Micromotor Handpiece
- Foot Control • Handpiece Stand

- Power Supply: AC 230 V 50/60 Hz
- Weight: 2.8 kg
- Dimensions: W 167 x D 247 x H 77 mm
- Max. Torque: 8.7 Ncm (Torque) / 6.0 Ncm (Compact)

MODEL **Torque Model: UMXL-GT**

ORDER CODE Y141469

MODEL **Compact Model: UMXL-GC**

ORDER CODE Y141470

RRP \$2932

SALE PRICE

\$2,199



NSK Ultimate XL-D Bench Top Tower

Contents

- Control Unit • Micromotor Handpiece
- Foot Control • Handpiece Stand

- Power Supply: AC 230 V 50/60 Hz
- Weight: 3.3 kg
- Dimensions: W 81 x D 270 x H 230 mm
- Max. Torque: 8.7 Ncm (Torque) / 6.0 Ncm (Compact)

MODEL **Torque Model: UMXL-DT**

ORDER CODE Y141504

MODEL **Compact Model: UMXL-DC**

ORDER CODE Y141505

RRP \$2932

SALE PRICE

\$2,199



NSK Ultimate XL-KT Knee Control

Contents

- Control Unit • Micromotor Handpiece
- Handpiece Stand

- Power Supply: AC 230 V 50/60 Hz
- Weight: 3.2 kg
- Dimensions: W 90 x D 270 x H 230 mm
- Max. Torque: 8.7 Ncm (Torque) / 6.0 Ncm (Compact)

MODEL **Torque Model: UMXL-KT**

ORDER CODE Y141495

MODEL **Compact Model: UMXL-KC**

ORDER CODE Y141496

RRP \$2932

SALE PRICE

\$2,199



NSK Ultimate XL-F Floor Type

Contents

- Control Unit • Micromotor Handpiece
- Handpiece Stand

- Power Supply: AC 230 V 50/60 Hz
- Weight: 3.8 kg
- Dimensions: W 230 x D 292 x H 113 mm
- Max. Torque: 8.7 Ncm (Torque) / 6.0 Ncm (Compact)

MODEL **Torque Model: UMXL-FT**

ORDER CODE Y141513

MODEL **Compact Model: UMXL-FC**

ORDER CODE Y141514

RRP \$2932

SALE PRICE

\$2,199

Additional Motor Cables & Handpieces

| MODEL | ORDER CODE | SALE PRICE |
|-----------------------------------------|------------|------------------------------|
| UMXL-T Torque Motor Handpiece 1.2m Cord | Y141440 | \$1,291 RRP \$1435 |
| UMXL-T Torque Motor Handpiece 2.0m Cord | Y141441 | \$1,873 RRP \$2082 |
| UMXL-C Compact Motor Handpiece 1.2m | Y141444 | \$1,291 RRP \$1435 |
| UMXL-C Compact Motor Handpiece 2.0m | Y141445 | \$1,873 RRP \$2082 |
| NSK Presto Aqua II Handpiece Only | T816 | \$1,487 RRP \$1652 |



Presto Aqua II

The PRESTO AQUA II offers adjustable water coolant options, ranging from fine mist to higher volume sprays, ideal for various materials. It features easy water bottle removal, can run independently or connect to mains water, and includes a dust-proof handpiece with a swivel, twist-free design.

MODEL **PRESTO AQUA II**

ORDER CODE Y150023

RRP \$3736

SALE PRICE

\$2,802



Contents

- PRESTO AQUA II Unit • PRESTO AQUA Handpiece • QD-J B3 Coupling • Hose for PRESTO AQUA II • Foot Control (AFC-40)
- Secondary Hose for PRESTO

- Water Coolant Spray
- Swivel Twist-free Handpiece
- Built-in Dust-proof System
- Simple & Easy Removal of Water Bottle