



Special Offers

JANUARY-MARCH 2025

Turbines

Experience Power

FREE TRIAL

Scan/click here to get started



Ti-Max Z

HIGH PERFORMANCE TURBINES

2 + 1 FREE*

NEW ULTIMATE POWER for quicker preps and cutting Zirconia

- Lightweight titanium body with Duragrip coating
- Also available in fittings for KaVo®, Dentsply Sirona®, W&H® (Z990/Z890 priced at RRP \$2,074, Z micro priced at RRP \$1,591)



MODEL **Z990L**
ORDER CODE **PA23730001**

- Optic
- **ULTIMATE Power** 30-44w**
- Std head: Ø 12.1 x 12.7 mm

**At 3 - 4.2 bar

RRP **\$1,981**



MODEL **Z890L**
ORDER CODE **PA23720001**

- Optic
- Power 28-41w**
- Std head: Ø 10.6 x 10.8 mm

RRP **\$1,981**



MODEL **Z micro (ultra-mini head)**
ORDER CODE **PA2357**

- Optic
- Power 20w
- Micro head: Ø 9 x 10.8 mm
- For Short Shank burs

RRP **\$1,381**

*Free item of equal or lesser value.

S-Max M

3 + 1 FREE*

- Stainless steel
- Enhanced power
- Quattro water spray



MODEL **M900L**
ORDER CODE **P1254**

- Optic
- Power 26w
- Std head: Ø 12.1 x 13.3 mm

RRP **\$1,381**

- DIY cartridge change
- Clean head system



MODEL **M800L**
ORDER CODE **P1253**

- Optic
- Power 23w
- Mini head: Ø 10.6 x 12.4mm

RRP **\$1,381**

- Also available in fittings for KaVo®, Dentsply Sirona®, W&H® (priced at RRP \$1,591)



MODEL **Pico**
ORDER CODE **P1140**

- Optic
- Power 9w
- Ultra-mini head: Ø 8.6 x 9mm
- Used only with short shank burs

RRP **\$1,381**

*Free item of equal or lesser value.

All prices as listed. Valid until 31 March 2025. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.



Ti-Max Z

3 + 1 FREE*

- Lightweight external titanium body
- **Optic with internal water**
- Small head/narrow neck
- **Unique spray switch on Z95L, Z45L**



MODEL **Z95L**
ORDER CODE **C1034SW**

- Optic
- 1:5 Increasing+Spray switch
- Head size Ø 8.9 x 13.4mm

RRP \$2,121



MODEL **nano95LS**
ORDER CODE **C1099**

- Optic
- 1:5 Increasing
- 10% shorter & lighter than Z95L
- Head size Ø 8.9 x 13.4mm

RRP \$1,909



MODEL **Z45L**
ORDER CODE **C1064**

- Optic
- 1:4.2 Increasing+Spray switch
- 45° head
- Head size Ø 9.6 x 14.5mm
- For regular or extra long (25mm) FG burs

RRP \$2,263



MODEL **Z25L**
ORDER CODE **C1038**

- Optic
- 1:1 Direct Drive
- Head Size Ø 8.9 x H 13.7 mm
- Single Water Spray
- For CA burs

RRP \$1,491

*Free item of equal or lesser value.

S-Max M

3 + 1 FREE*

- Stainless steel
- **All with internal water**
- Choice of optic & non-optic



MODEL **M95L**
ORDER CODE **C1023**

- Optic
- 1:5 Increasing
- Head size 10 x H.14.6mm

RRP \$1,763



MODEL **M25L**
ORDER CODE **C1024**

- Optic
- 1:1 Direct Drive

RRP \$1,353



MODEL **M15L**
ORDER CODE **C1025**

- Optic
- 4:1 Reducing

RRP \$1,647



MODEL **M65**
ORDER CODE **H1008**

- Non-Optic
- 1:1 Straight

RRP \$818

MODEL **M95**
ORDER CODE **C1026**

- Non-Optic
- 1:5 Increasing
- Head size 10 x H.14.6mm

RRP \$1,405

MODEL **M25**
ORDER CODE **C1027**

- Non-Optic
- 1:1 Direct Drive

RRP \$949

MODEL **M15**
ORDER CODE **C1028**

- Non-Optic
- 4:1 Reducing

RRP \$1,334

*Free item of equal or lesser value.

Contra Angles

SCAN/CLICK
for more info



FXm

4 + 1 FREE*



MODEL **FX65m**
ORDER CODE **H1014001**

- Non-Optic Straight
- 1:1 Direct Drive
- Push button chuck
- No internal water

RRP | \$461



MODEL **FX25m**
ORDER CODE **C1052001**

- Non-Optic Contra-angle
- 1:1 Direct Drive
- Push button chuck
- No internal water

RRP | \$462



MODEL **FX57m**
ORDER CODE **C1054001**

- Non-Optic Prophy Contra-angle
- 4:1 speed reducing
- Max speed 1,000 rpm
- Brushes/Cups: Screw in

RRP | \$508



MODEL **FX22**
ORDER CODE **C1129091**

- Non-Optic Contra-angle
- 1:1 Direct Drive
- Latch-type chuck
- No internal water

RRP | \$279

*Free item of equal or lesser value.

Maintenance

SCAN/CLICK
for more info



iCare HANDPIECE LUBRICATION UNIT

- Precise dosing – less waste, better care
- Connects to 6mm airline, 5-6 bar

MODEL **iCare C2 Type**
ORDER CODE **Y1002787**

RRP | \$5,489

MODEL **iCare C3 Type**
ORDER CODE **Y1002788**

RRP | \$5,863



Buy iCare C2 or C3 and get **\$2,304** worth of **FREE Handpieces + adaptors**

Handpieces include: 1 x S-Max High-Speed, 1 x FX65M and 1 x FX25M, PLUS Free High-Speed Adaptors: 2 x Adaptors for iCare C2, 1 x Adaptor for iCare C3.



Surgic Pro and Surgic Pro2

MICROMOTORS FOR IMPLANTS / ORAL SURGERY

Surgic Pro Includes:

- Control unit
- Motor/Cord
- 20:1 Implant contra angle
- Foot control
- 2 year warranty



MODEL **Surgic Pro (optic)**
ORDER CODE **Y1001933**

RRP | \$8,726

Surgic Pro2 Includes:

- Control unit
- Motor/Cord
- 20:1 Implant contra angle
- Wireless foot control
- Bluetooth connectivity to VarioSurg & Osseo100+
- 2 year warranty



MODEL **Surgic Pro2 (optic)**
ORDER CODE **Y1500511**

RRP | \$8,726

Buy Surgic Pro or Surgic Pro 2 Optic and choose from
FREE Additional Micromotor & Cable OR ISD900 OR 2 x Ti-Max X-SG65L

Micromotor



Optic Micromotor Cable for
Surgic Pro2 Optic

MODEL **SGL80M**
ORDER CODE **Y1004211**

RRP | \$4,464

Optic Micromotor Cable for
Surgic Pro Optic

MODEL **SGL70M**
ORDER CODE **E1023001**

RRP | \$4,464

iSD900



MODEL **iSD900**
ORDER CODE **Y1001358**

RRP | \$3,604

Ti-Max X-SG65L



MODEL **Ti-Max X-SG65L**
ORDER CODE **H1009**

RRP | \$2,371



VarioSurg3

VARIOSURG 3 ULTRASONIC BONE SURGERY SYSTEM

50% OFF for existing SurgicPro or SurgicPro2 users

Features:

- Control unit
- LED handpiece with 2m cord
- FC-78 foot control
- Sterilisation cassette
- Handpiece stand
- Irrigation tube (5 pieces)
- Other accessories

MODEL **VarioSurg 3**
(without foot control)
ORDER CODE **Y1002248**

RRP **\$10,060**

MODEL **VarioSurg 3**
ORDER CODE **Y1002726**

RRP **\$12,468**



iCart

15% OFF RRP

The NSK iCart Duo holds both the VarioSurg3 & Surgic Pro2. Use one foot control and link together with the Link Module.

MODEL **iCart**
ORDER CODE **S9090**

RRP **\$2,354**

15% OFF **\$2,001**

MODEL **Link Module (For Surgic Pro 2)**
ORDER CODE **Z1402001**

RRP **\$370**



VarioSurg Tips

2 + 1 FREE*

SCAN/CLICK
to check out our full
range of VarioSurg tips



Tip Holder sold separately

*Free item of equal or lesser value.

Oral Hygiene

Prophy-Mate neo

SUPRAGINGIVAL POWDER THERAPY

ORDER CODE **Y135029** for NSK PTL coupling

- Compact, lightweight & portable
- Includes 60° and 80° nozzles
- 360° rotation
- Includes 6 sachets of FLASH pearl® powder
- Also available in fitting for KaVo®, Dentsply Sirona®, W&H® (priced at RRP \$1,809)

RRP **\$1,696**



SCAN/CLICK
for more info



Buy both the Prophy-Mate and the Perio-Mate and get 25% off RRP for each



+ **FREE** box of **FLASH pearl** or **SOFT pearl Powder**



FLASH pearl

- Calcium carbonate (54µm) based powder for stain removal

4 x 300 g bottles

ORDER CODE **Y900698**

RRP **\$313**



SOFT pearl

- Glycine based powder (65µm), ideal for moderate staining. Water soluble, making clean up easier

4 x 250 g bottles

ORDER CODE **Y1500721**

RRP **\$326**

NEW



SCAN/CLICK
for more info

Perio-Mate

SUBGINGIVAL POWDER THERAPY / BIOFILM REMOVAL

CODE **Y1002654** for NSK coupling

- Subgingival version of ProphyMate
- Eradicates biofilm in just 20 seconds
- Powder flow adjuster for delicate pockets
- Includes flexible plastic nozzle tips for deeper pockets
- Packed in a handy case with 160g powder
- Also used supragingivally on moderate staining and around composite restorations/ orthodontic brackets
- Also available in fitting for KaVo®, Dentsply Sirona®, W&H®

RRP **\$1,696**

SCAN/CLICK
for more info



Buy both the Prophy-Mate and the Perio-Mate and get 25% off RRP for each



+ **FREE** box of **Perio-Mate Powder and Nozzle Tips**

FREE



MODEL Perio Mate Powder

ORDER CODE **Y900938**

- Glycine based powder (25µm)
- 2 x 160 g bottles

RRP **\$313**




MODEL Perio Mate Nozzle Tips

ORDER CODE **Y1002741**

- Pack of 40

RRP **\$125**



POWDER THERAPY GUIDE



SCAN/CLICK
to download.

Oral Hygiene

Varios Combi Pro IPIEZO ENGINE®

TWO ORAL CARE SYSTEMS IN ONE!

FREE Perio Set with Varios Combi Pro Basic Set

- Combines ultrasonic treatments and powder therapy in one unit
- Comes with a detachable Prophy chamber
- Detachable Perio chamber available as an optional extra

MODEL **Varios Combi Pro Basic Set**

ORDER CODE **Y1004194**

RRP \$7,856

MODEL **Varios Combi Pro Perio Set**

ORDER CODE **Y1003042**

RRP \$2,839

Ultrasonic

SCAN/CLICK
for more info



Perio & Prevention

SCAN/CLICK
to check out our full
range of Scaler tips




Scaler Tips

TIPS FOR VARIOS ULTRASONIC SCALER

2 + 1 FREE*

NSK offers a Wide Choice of Varios Ultrasonic Scaler Tips for Clinical Applications.

*Free item of equal or lesser value.





Endo-Mate

ENDO-MATE DT2

The ENDO-MATE DT2, compatible with major Ni-Ti file brands, features a pre-programming function for efficient use of various file protocols. A lightweight, compact control unit with a rechargeable battery for enhanced mobility. The slim, lightweight handpiece ensures comfortable handling and high performance.

Includes:

- Control unit
- Motor and cord
- MP-F20R head (20:1)
- Handpiece stand
- AC adaptor

Specifications:

- Power Supply 230V 50/60 Hz and Battery pack
- Number of Programs 9
- Max Torque 7 Ncm (when using 20:1 head)
- Speed range 100-13,000 min⁻¹ (With 20:1, 1:1 geared heads)
- Charge/Continuous use 5 h / 2 h (depending on workload)
- Dimensions W 93 x D 147 x H 100 (mm)
- Weight Unit 440 g / Motor handpiece and Cord 95g

MODEL **Endo-Mate DT2**

ORDER CODE **Y1004230**

RRP **\$2,348**



ENDO-MATE TC2

The Cordless ENDO-MATE TC2 Endodontic Handpiece enhances safety and efficiency with advanced torque control, auto-reverse, and a slim design for improved visibility. It features programmable settings, dual alarm sounds, and compatibility with most Ni-Ti files. An integrated LCD panel and the ability to combine with an Apex locator facilitate precise treatments. The device includes a quick 90-minute non-contact charger for fast recharging.

ENDO-MATE TC2 includes:

- MP-F16R Head
- ENDO-MATE TC2 Motor
- Battery Charger

ENDO-MATE TC2 MPA includes:

- MPA-F16R Head
- ENDO-MATE TC2 Motor
- Battery Charger

Specifications:

- Power Supply 230V 50/60 Hz and Battery pack
- Number of Programs 9
- Max Torque 7 Ncm (when using 20:1 head)
- Speed range 100-13,000 min⁻¹ (With 20:1, 1:1 geared heads)
- Charge/Continuous use 5 h / 2 h (depending on workload)
- Dimensions W 93 x D 147 x H 100 (mm)
- Weight Unit 440 g / Motor handpiece and Cord 95g



Buy ENDO-MATE DT2 or ENDO-MATE TC2 and get **iPexII 50% OFF**

iPexII

APEX LOCATOR

MODEL **iPexII**

ORDER CODE **Y1002208**

RRP **\$1,696**

SALE PRICE **\$848***



MODEL **ENDO-MATE TC2**

ORDER CODE **Y1001027**

RRP **\$2,263**

MODEL **ENDO-MATE TC2 with MPA**

ORDER CODE **Y1001028**

RRP **\$2,716**



* When you purchase either an ENDO-MATE TC2 or ENDO-MATE DT2 Complete set.

CONTACT NSK

Australia

Sales & Education Manager
Mariam Dalzotto
mariam@nskoceania.com.au
0427 833 684

Queensland

Business Development Manager
Sam Vaizie
sam@nskoceania.com.au
0439 260 006

NSW North

Business Development Manager
Leah Nicholas
leah@nskoceania.com.au
0419 251 216

NSW South & ACT

Business Development Manager
Janette Carroll
janette@nskoceania.com.au
0439 199 979

Victoria & Tasmania

Business Development Manager
Dr Olivier Nguyen
olivier@nskoceania.com.au
0456 551 087

West & South Australia

Business Development Manager
Amar Mandalia
amar@nskoceania.com.au
0459 908 999